

# MATSUNAGA

Wheelchair, Walker, Walking Cane, Stretcher, Bathing Aid

FULL LINE CATALOGUE vol.1

アサツの

支えになりたい



Matsunaga Manufactory Co., Ltd.

## ISO9001 Certified

The Matsunaga brand is ISO9001 certified to guarantee high reliability in terms of the quality, shipping, and cost of our products.

## Matsunaga's New Manufacturing Office Shanghai Factory

We are ISO9001/ISO13485 certified to ensure high quality manufacturing practices. We promise higher quality products than ever before, endeavouring to deliver higher functionality and greater value to our customers.

## JIS Mark Certified for Assistive Products

### JIS Mark Certified Models

DM-81/AR-101/AR-111/AR-201B/AR-211B/ AR-301/AR-311/  
AR-501/AR-601/\*AR-901(Self-Propelled, Attendant-Propelled)/  
MW-SL11B/MW-SL21B/MW-SL31B/MW-SL41B/NEXT-11B/  
NEXT-21B/NEXT-31B/NEXT-41B/\*NEXT-51B/\*NEXT-61B/  
NEXT-10B/NEXT-12B/NEXT-20B/NEXT-22B/\*NEXT-50B/  
\*NEXT-60B/GRC-11B/GRC-21B/GRC-31B/GRC-41B

Standard Number JIS T9201

\* The JIS Mark is given at the time of manufacturer shipment.

Certain models may not be JIS compliant with optional parts installed.



そうそう、  
マツナガセイサクショ  
50年前、障害をもった方に  
1台1台手作りの車いすを作り始めたのがキッカケだったんだっけ。  
下を向いて過ごしていた方が上を向き、  
止まっていた地点から前へ、一步一步少しずつ一緒に。  
セカンドチョイスで構わない、  
キミの傍にずっとマツナガセイサクショ。



# Table of Feature Icons

Please keep the catalogue flipped open to this page as a reference when going through each category of wheelchairs.



Please keep this page flipped open.

	<b>Fixed Backrest</b>	This fixed type of backrest frame is highly durable.
	<b>Foldable Backrest</b>	The backrest can be folded. This features makes it easier to load the wheelchair onto cars.
	<b>Forward-Folding Backrest</b>	The backrest can be folded forward. This features makes it easier to load the wheelchair onto cars.
	<b>Detachable Armrests</b>	The armrest pipes are removable. This feature is convenient when transferring the user onto beds.
	<b>Vertically Adjustable Armrests</b>	The armrests can be slid up and down to an optimal position to support the user's elbows.
	<b>Flip-Up Armrests</b>	The armrests can be easily flipped up with one touch.
	<b>Flip-Up &amp; Detachable Armrests</b>	The armrests can be flipped up with one touch. The armrest pipes are removable.
	<b>Height Adjustable Armrests</b>	The armrest height can be adjusted to an optimal position according to the user's physique.
	<b>Armrest Invertible Pads</b>	The direction of the armrest pads can be inverted from left to right, and vice versa. The armrests can be extended forward to support the user when standing up.
	<b>Adjustable Backrest Upholstery</b>	This feature allows long durations of wheelchair use without experiencing fatigue. The backrest upholstery belts can be adjusted according to the shape of the user's back.
	<b>Adjustable Backrest · Seat Upholstery</b>	This feature allows long durations of wheelchair use without experiencing fatigue. Each backrest and seat upholstery belts can be independently adjusted according to the user's physique.
	<b>Swing-Away</b>	The footrest · legrest units can be swung in and out, or detached with one touch. This is convenient when positioning the wheelchair closer to beds or when foot propelling.
	<b>Swing-In-and-Away</b>	The footrest · legrest units can be swung in and out, or even straight up for detachment. This feature is convenient for use in narrow spaces.
	<b>Elevating</b>	The angle of the footrest · legrest units can be adjusted with one touch. This feature is suitable for users with foot rigidity or deformities.
	<b>Elevating Separate</b>	The elevating and reclining features of the wheelchair function separately.
	<b>Elevating Interlocked</b>	The elevating and reclining features of the wheelchair are interlocked.
	<b>Reclining (Manual)</b>	The angle of the backrest can be adjusted using a lever.
	<b>Full Reclining</b>	The reclining angle can be adjusted to an almost fully flat mode. This feature is suitable for users on bed rest. The wheelchair can also be used as a simple stretcher.
	<b>Tilt</b>	The seat angle can be adjusted to alleviate strain on the user in sitting position. This also prevents the user from sliding forward.
	<b>Tilt &amp; Reclining</b>	The user's body weight can be more evenly dispersed by adjusting both the seat and backrest angles.
	<b>OASIS Double Reclining</b>	This allows the user to maintain a good posture more easily. Reclining only engages the pelvic area to tilt, thereby preventing the user from sliding forward.
	<b>Drumbrakes</b>	The brake can be easily operated by the attendant. It is even effective on slopes and streets wet from the rain.
	<b>Anti-Tip Bars (Fall Prevention)</b>	The metal fixture is designed for fall prevention to keep the user's center of gravity from shifting backwards when going uphill.
	<b>Footbrakes</b>	The foot-operated brakes can be activated with one touch, enabling braking without taking either hands off the wheelchair.
	<b>Duralock</b>	This brake features drumbrakes coupled with footbrakes to enable stable braking performance.
	<b>Antibacterial</b>	The frame pipes have antibacterial properties for clean use at all times.
	<b>3D Frame</b>	We trimmed down the total width of the wheelchair while keeping the seat width the same. This feature is suitable for home use.
	<b>Super Casters</b>	These non-puncture tires are made of polyurethane foam with high shock absorption performance. They reduce vibration on even terrain or bumpy roads. This feature is suitable for outdoor use.
	<b>Puncture Proof Tires</b>	There is no need to pump air into the tires. This feature is suitable for users unable to frequently inflate tires.
	<b>Safety Belt Included</b>	This feature prevents falling off the wheelchair and bad posture.
	<b>Detachable Headrest</b>	The headrest can be detached when assisting from behind to alleviate strain on the attendant.
	<b>Height Adjustable Headrest</b>	The height of the headrest can be adjusted according to the user's head position.
	<b>Height Adjustable Seat (Manual)</b>	The seat can be slid up and down to an optimum position.
	<b>Position &amp; Angle Adjustable Footrests</b>	The position and angle of the foot plates can be adjusted accordingly. This is suitable even for users suffering from malformation of feet joints (equinus contracture).

# How to Navigate This Catalogue

Please make sure to read this section first to navigate this catalogue properly.

**1 Table of Feature Icons**

Please keep the catalogue flipped open to this page as a reference when going through each category of wheelchairs.

**2 Lightweight Multi-Functional Wheelchair**

**3 Lightweight Flip-Up Armrest & Swing-Away Footrests**

**4 Front Seat Height 430 mm**

**5 MW-SL31B**

**6**

**7 Lightweight Flip-Up Armrest & Swing-Away Footrests**

**8 MW-SL41B**

**9**

**10**

## 1 Table of Feature Icons

Please keep the List of Feature Icons flipped open while going through this catalogue. Here you can easily refer to descriptions of the features mentioned on each product page.

## 2 Category

The products are categorized according to series and model here.

## 3 Model

The product type is indicated here.

## 4 Front Seat Height (Front Seat to Floor Height)

This is the distance between the floor and the wheelchair seat (top of the seat frame).

## 5 Product Photo

The product type can be confirmed in the photo.

## 6 Seat Upholstery Colors

The range of seat upholstery colors available for the product is shown here.

## 7 Features

The icons specify the features of each product.

## 8 Manufacture's Suggested Retail Price

The price and the taxability status of the product is indicated here.

## 9 Model Number

This number is used to classify the product.

## 10 Size Chart

The dimensions of each product part is indicated here. Please refer to this chart along with the model number.

## Choosing the Size

To choose a suitable wheelchair for yourself...

First and foremost, the wheelchair size must be right for you. It's strenuous for the user if the wheelchair is too big or too small. Please check out the tips below.

### Front Seat Height (Front Seat to Floor Height)

As a rule of thumb, a comfortable seat height is 5 cm higher than the length of the user's lower leg. However, this should be only 1~2 cm higher for users who foot propel. Also, users with limited knee mobility should opt for higher seat heights. (Please take the thickness of the seat cushion into consideration.)

### Seat Width

(Distance Between Armrest Pipes)  
Please allow adequate clearance for one hand to fit on either side of the user's buttocks.

### Handgrip Height

The handgrips should be positioned at a height easiest for the attendant to push with elbows slightly bent (generally the height of the attendant's belly button).

### Armrest Height


As a rule of thumb, this should be 2~3 cm higher than the height of the user's elbows bent at 90° relative to the seat. Select a height neither too high nor too low. (Please take the thickness of the seat cushion into consideration.)

### Footrest Height

The footrests should be more than 5 cm from the floor, with the user's thighs slightly touching the seat. The distance between the seat and thighs should be neither too far apart nor too close together.

### Backrest

For those able to self-propel, the backrest should come up to the bottom of the user's shoulder blades.



Seat Width
Hip Width
+
4~5 cm of Clearance

Backrest Height
Distance from Seat to Bottom of Scapula
+
4~5 cm of Clearance

# C O N T E N T S

## New Lineup

Grace Core Standard — 6      Grace Core - Multi-Functional — 7

5

New Lineup

## NEXT CORE Series

NEXT CORE - Adjust	12
NEXT CORE - Minimo	13
NEXT CORE - Wide	14
NEXT CORE - Multi-Functional	15
NEXT CORE - Standard	16
NEXT CORE - Grand • Petit	17
NEXT CORE - Spin	18

Choosing NEXT CORE Based on Height & Features	19
3 Tips for Sitting Stability on NEXT CORE	20
How to Adjust the Backrest Upholstery (Recommendation)	21
NEXT CORE Series Option Chart	22

11

NEXT CORE  
Series

## Standard • Multi-Functional Wheelchairs

DM Series	24	AR Series	28	Aerial Series	38
CM Series	26	MY Series	36		

23

Standard •  
Multi-Functional

## Reclining • Tilt Wheelchairs

My Tilt Series	44	OASIS Reclining Series	50	CM Series	54	Escort	56
OASIS Positive Series	51			MW Series	55	Elise	57

43

Reclining • Tilt

## Child Series

LUUK	60	Mignon	63	Tight Turn Jr.	64
MB Series	62	Leveret	64		

59

Child

## Special Features Series

Auto Break Series	66	One-Arm Drive Wheelchairs	68
MRI Series	67		

65

Special  
Features

## Nursing Care Products

Compact Wheelchairs • Walkers • Walking Canes • Bathing Aids, etc.

Outdoor Walker	70	Indoor Walker	72	Walking Cane	74
Compact Wheelchair	71	Walker	73	Bathing Aids	75

69

Nursing Care  
Products

## Wheeled Stretcher • Stretcher

Aluminum Lightweight Wheeled Stretcher	80	Stretcher	85
MRI-Compatible Wheeled Stretcher	84		

79

Wheeled Stretchers  
• Stretchers

Product Lineup — 2

Index — 88

Company Profile — 90



## Product Lineup

Matsunaga aims to become a comprehensive manufacturer of medical·assistive products.

### Manual Wheelchairs

#### New Lineup



Grace Core  
▶P6~P9



NEXT CORE Series  
▶P12~P22



DM Series  
▶P24~P25



CM Series  
▶P26~P27



AR Series  
▶P28~P35



MY Series  
▶P36

#### Standard·Multi-Functional

#### Ultra Lightweight



Aerial Series  
▶P38~P40

#### Tilt & Reclining



My Tilt Series  
▶P44~P49



OASIS Positive &  
Reclining Series  
▶P50~P53



CM Series  
▶P54



MW Series  
▶P55



Escort  
▶P56



Elise  
▶P57

#### Child Series



LUUK  
▶P60~P61



MB Series  
▶P62



Mignon  
▶P63



Leveret  
▶P64



Tight Turn Jr.  
▶P64

### Compact Wheelchairs · Walkers · Walking Canes · Bathing Aids

#### Outdoor Walker



OLive  
▶P70

#### Compact Wheelchair



MV Series  
▶P71

#### Indoor Walker



Hop Step Series  
▶P72

#### Reciprocal·Standard Walking Frame



CMS Series  
▶P73

#### Walking Cane



Quad Cane Lofstrand  
Crutch  
▶P74

#### Bathing Aids



Coral Series  
▶P75~P77

#### Special Features Series



Auto Break Series  
▶P66



MRI Series  
▶P67



One-Arm Drive Wheelchairs  
▶P68

#### Semi-Module Model



NEXT-51B/61B  
Self-Propelled  
Attendant-Propelled  
▶P12



NEXT-50B/60B  
Self-Propelled  
Attendant-Propelled  
▶P13



AR-901  
Self-Propelled  
Attendant-Propelled  
▶P33~P34



AR-911S  
Self-Propelled  
Attendant-Propelled  
▶P35

### Hospital · Care Facility Equipment


#### Aluminum Lightweight Wheeled Stretcher



AL Series  
▶P80~P83

### Emergency Equipment

#### Stretchers



2-Fold·4-Fold Stretchers with Legs  
▶P85



Assist Stretchers  
▶P86



# Bringing Matsunaga's Hospitality to the World.

From airport counters to onboard aircrafts, our barrier-free wheelchair morph made of 100% resin is the way to go. It marks a new beginning for Matsunaga, from a Japan-centered Matsunaga to a global Matsunaga. We'll always be there for each and every person.

Product information  
available on YouTube



**m o r p h**

MATSUNAGA MANUFACTORY CO., LTD.



**GOOD DESIGN AWARD 2016  
GOLD AWARD**

(PIIF) Onogi Science and Technology  
Promotion Foundation  
Winner of the Top Silver Award in 2016



**GERMAN  
DESIGN  
AWARD  
WINNER  
2018**



Winner of the  
German Design Award

# new lineup



New Lineup



Grace Core Standard



GRC-11B (Self-Propelled)  
GRC-21B (Attendant-Propelled)



The long-awaited high-grade model of  
our NEXT CORE Series is finally available!



\* May not be JIS  
compatible with  
optional parts  
installed.

Load Capacity  
100kg



Seat Upholstery Colors



M-11 M-12 M-13

Grace Core Standard Model • Self-Propelled

Front Seat Height **430 mm**

\* Seat cushion not included.

**GRC-11B** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥125,000**




Grace Core Standard Model • Attendant-Propelled

Front Seat Height **430 mm**

\* Seat cushion not included.

**GRC-21B** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥125,000**



Grace Core Standard Model • Attendant-Propelled with Duralocks (Foot Brakes)

Front Seat Height **430 mm**

\* Seat cushion not included.

**GRC-21BD** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥145,000**



Duralock Design  
(Foot Brakes)

Unit: mm

Product Number	Caster x Rear Wheel	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Length	Total Height	Total Width	Weight
<b>GRC-11B</b>	5×22	430	400/420	*380	415	250	560	290~350	950	930	565	13.8kg
<b>GRC-21B</b>	5×16	430	400/420	*380	415	250	560	290~350	950	930	565	12.4kg
<b>GRC-21BD</b>	5×16	430	400/420	*380	415	250	560	290~350	950	930	565	13.7kg

\* The seat depth can be adjusted to 380mm or 340mm.(The factory setting of the seat depth is 380mm.)

Grace Core - Multi-Functional



GRC-31B (Self-Propelled)  
GRC-41B (Attendant-Propelled)



The Multi-Functional Model of Our Grace Core Series  
with Flip-Up • Swing-In • Swing-Away Features!



\* May not be JIS  
compatible with  
optional parts  
installed.

Load Capacity  
100kg



Seat Upholstery Colors




M-11 M-12 M-13

Grace Core Multi-Functional Model • Self-Propelled

Front Seat Height **430 mm**

\* Seat cushion not included.

**GRC-31B** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥148,000**




Grace Core Multi-Functional Model • Attendant-Propelled

Front Seat Height **430 mm**

\* Seat cushion not included.

**GRC-41B** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥148,000**



Grace Core Multi-Functional Model • Attendant-Propelled with Duralocks (Foot Brakes)

Front Seat Height **430 mm**

\* Seat cushion not included.

**GRC-41BD** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥168,000**



Duralock Design  
(Foot Brakes)

Unit: mm

Product Number	Caster x Rear Wheel	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Length	Total Height	Total Width	Weight
<b>GRC-31B</b>	5×22	430	400/420	*380	415	250	560	290~350	950	930	565	15.5kg
<b>GRC-41B</b>	5×16	430	400/420	*380	415	250	560	290~350	950	930	560	14.1kg
<b>GRC-41BD</b>	5×16	430	400/420	*380	415	250	560	290~350	950	930	565	15.4kg

\* The seat depth can be adjusted to 380mm or 340mm.(The factory setting of the seat depth is 380mm.)



## Our Airless Hybrid Tires!

Our hybrid tires have inner tubes made of polyurethane foam with excellent shock absorbercy instead of conventional air-filled tubes. This standard feature makes pumping air into the tires unnecessary. Easy to operate, the brakes last longer to ensure user safety.

Shrink-resistant high density polyurethane foam tube



Here's a detailed video of the comprehensive standard features of Grace Core!


[What a Difference! Our Standard Wheelchair Features!](#)

## Our Easy Dry Cushion with Excellent Breathability!


The seat upholstery is made of our mesh fabric specially designed to resist pilling. The cushion is made of foam. Highly breathable and moisture-wicking, the seat is comfortable to sit in.



Specially Designed Mesh Fabric



Foam Cushion



Easy dry cushion  
Weight: 0.5kg

## Our New Backrest Mechanism for Comfort Even After Hours of Sitting!


### ◆ Our Standard 56cm Backrest is the Tallest in the Industry!

Grace Core Backrest Height 56cm

NEXT CORE Backrest Height 46cm




The backrest is 10 cm higher than in our NEXT CORE series. The wider backrest supports the back to offer greater head stability.



Grace Core's Back Frame

NEXT CORE's Back Frame

The backrest frame curves further towards the back than our NEXT CORE series, making it more relaxing to lean back in.



Here's a detailed video of the secrets of Grace Core backrests!

[What a Difference! The Secrets of our Backrests!](#)


### ◆ Our Backrest Upholstery with Triangular Pad for Outstanding Back Support & Maneuverability!



Gaps form around the shoulders when applying pressure to the back frame. The triangular pad attached to the upper section of the backrest cover seat supports the back, particularly the spinal column, to maintain the gap.




Part with Triangular Pad for Support



## Comfortable even after hours of sitting!

The high backrest fits the entire back for proper support. This dramatically reduces the pressure on the back and seat areas.




[Grace Core]

[Conventional Wheelchair]

Users can move their arms more easily thanks to the triangular pad!

The backrest pushes on the shoulder blades. Users cannot pull back their arms...

The backrest with triangular padding properly supports the spinal column. The backrest does not push on the shoulder blades, so users can easily move their arms!




[Grace Core]

[Conventional Wheelchair]

## Further Upgrades to Grace Core's Adjustable Backrest Upholstery!

### ◆ Better backrest adjustment for a better fit

In our NEXT CORE Series, there is a gap between the 2nd and 3rd belts from the top. This requires careful placement of the belts.



Backrest Belts of NEXT CORE



Backrest Belts of Grace Core

Backrest belts with less gaps in between and a frame which better fits the curvature of backs **makes it possible to adjust the backrest like a pro!**

## Comes newly equipped with innovative gluteal support belts!

### ◆ The gluteal support belts prevents gaps from forming

It is possible to minimize pressure on the hipbone by reducing the gap between the posterior pelvis and the seat to support the entire buttocks. This also prevents the pelvis from tilting backwards.





Image of gluteal support belts

Pelvic support  
Gluteal support

Just like NEXT Core, Grace Core comes with pelvic support to prop up the pelvis from all directions. Our new feature - the gluteal support belts - similarly provides the lower pelvis further support.




Here's a detailed video of the secrets of Grace Core backrests!

[What a Difference! Tips on Adjusting Backrests!](#)

## Easily Adjust the Seat Depth on the Spot!

The seat depth can be easily adjusted to 38cm or 34cm without any tools. Adjust the seat depth to a size according to the user's height and usage right on the spot.



Seat Depth 38cm

Seat Depth 34cm



Here's a detailed video of the secrets of Grace Core seats!

[What a Difference! The Secret of our Seats!](#)



[Grace Core]



[Conventional Wheelchair]



Here's a detailed video of the secrets of Grace Core footrests!

[What a Difference! The Secret of our Footrests!](#)

## Our 90 Degree Footrests!

It is easy to sit in the wheelchair with the footrests positioned closer in towards the user like in a chair.

## GRC Option Chart

	Walking Cane Holder (NEXT)	Posterior IV Pole (Standard)	Oxygen Tank Holder (NEXT)	Foldable Oxygen Tank Holder (NEXT)	Anti-tip Bars (Fall Prevention) (NEXT)	Extension Bar (5cm/10cm) (Standard)	Spoke Cover (NEXT)
GRC-11B	●	●	●	●	●	●	●
GRC-21B	●	●	●	●	●	●	×
GRC-31B	●	●	●	●	●	●	●
GRC-41B	●	●	●	●	●	●	×
GRC-21BD	×	×	×	×	×	●	×
GRC-41BD	×	×	×	×	×	●	×

\* Further discussion needed when combining several options

\* Please refer to P.22 for images of the optional parts when installed

## NOTE

---



## NEXT CORE Series






NEXT-51B HB/51B (Self-Propelled)  
NEXT-61B HB/61B (Attendant-Propelled)



A modular wheelchair supporting good posture  
with adjustable armrests • seat • footrests



Seat Upholstery Colors



NEXT-51B (Self-Propelled)

Height Adjustable Armrest



Convenient for use with  
commercially sold cushions!  
\* 25cm~31cm

Position & Angle  
Adjustable Footrests



Adjust to fit the user's  
height and physical  
condition!

Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled Front Seat Height 410 mm / 430 mm / 450 mm \* Seat cushion not included.  
NEXT-51B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥162,000 Caster x Rear Wheel [Inches] 6 x 22 (22 x 1 1/4) Weight 15.2kg

Attendant-Propelled Front Seat Height 410 mm / 430 mm / 450 mm \* Seat cushion not included.  
NEXT-61B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥162,000 Caster x Rear Wheel [Inches] 6 x 16 (16 x 1 3/8) Weight 14.2kg

Duralock Design (Foot Brakes)

Attendant-Propelled Front Seat Height 410 mm / 430 mm / 450 mm \* Seat cushion not included.  
NEXT-61B HBD Manufacturer's Suggested Retail Price (Non-Taxable) ¥182,000 Caster x Rear Wheel [Inches] 6 x 16 (16 x 1 3/8) Weight 15.5kg

Air Tire Design

Self-Propelled Front Seat Height 410 mm / 430 mm / 450 mm \* Seat cushion not included.  
NEXT-51B Manufacturer's Suggested Retail Price (Non-Taxable) ¥156,000 Caster x Rear Wheel [Inches] 6 x 22 (22 x 1 1/4) Weight 14.6kg

Attendant-Propelled Front Seat Height 410 mm / 430 mm / 450 mm \* Seat cushion not included.  
NEXT-61B Manufacturer's Suggested Retail Price (Non-Taxable) ¥156,000 Caster x Rear Wheel [Inches] 6 x 16 (16 x 1 3/8) Weight 13.6kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-51B HB	410/430/450	380/400/420	380	380/400/420	250~310	460	325~415	880/900/920	960	560
NEXT-51B										550
NEXT-61B HB	410/430/450	380/400/420	380	380/400/420	250~310	460	325~415	880/900/920	955	560
NEXT-61B										550


\* Footrest length: distance between the footrests and seat

\* The factory setting of the front seat height is 430mm.

\* Additional fees are required for adjustments to the front seat to floor height before shipping.

© Specifications are subject to change without prior notice for continual product improvement.

© The dimensions are based on the seat width set at 400mm.




NEXT-61B  
(Attendant-Propelled)

NEXT-50B HB/50B (Self-Propelled)  
NEXT-60B HB/60B (Attendant-Propelled)




An extra-low to low-seat adjustable modular  
wheelchair for use with cushions and foot propelling



NEXT-50B (Self-Propelled)

Seat Upholstery Colors



Height Adjustable Armrest



Convenient for use with  
commercially sold cushions!  
\* 22.5cm~28.5cm

Height Adjustable Seat



Adjust to fit the user!  
\* Front Seat Height 39cm/37cm/35cm

Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled Front Seat Height 390 mm / 370 mm / \*350 mm \* Seat cushion not included.  
NEXT-50B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥167,000 Caster x Rear Wheel [Inches] 6 x 20 (20 x 1 1/4) Weight 14.9kg

Attendant-Propelled Front Seat Height 390 mm / 370 mm / \*350 mm \* Seat cushion not included.  
NEXT-60B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥167,000 Caster x Rear Wheel [Inches] 6 x 16 (16 x 1 3/8) Weight 14.2kg

Air Tire Design

Self-Propelled Front Seat Height 390 mm / 370 mm / \*350 mm \* Seat cushion not included.  
NEXT-50B Manufacturer's Suggested Retail Price (Non-Taxable) ¥161,000 Caster x Rear Wheel [Inches] 6 x 20 (20 x 1 1/4) Weight 14.3kg

Attendant-Propelled Front Seat Height 390 mm / 370 mm / \*350 mm \* Seat cushion not included.  
NEXT-60B Manufacturer's Suggested Retail Price (Non-Taxable) ¥161,000 Caster x Rear Wheel [Inches] 6 x 16 (16 x 1 3/8) Weight 13.6kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-50B HB	390/370/350	400	360	365 / 345 / 325	225~285	460	285~375	880 / 860 / 840	940	560
NEXT-50B								*Changed to 5-inch casters.		
NEXT-60B HB	390/370/350	400	360	365 / 345 / 325	225~285	460	285~375	880 / 860 / *840	940	550
NEXT-60B								*Changed to 5-inch casters.		

\* Footrest length: distance between the footrests and seat

\* The factory setting of the front seat height is 390mm.

© Specifications are subject to change without prior notice for continual product improvement.

\* Refer to the chart of options available on P22.



NEXT-60B  
(Attendant-Propelled)

Chart of Wheelchair Sizes  
(NEXT-50B • 60B  
(Self-Propelled)(Attendant-Propelled))

Unit: mm

	Front Seat Height	Caster Size	Wheel Size
Self-Propelled	390	6-inches	20x11/4
	370	6-inches	20x11/4
	350	5-inches	20x11/4
Attendant-Propelled	390	6-inches	16x13/8
	370	6-inches	16x13/8
	350	5-inches	16x13/8

\* The 5-inch casters are solid cushioned casters.



A New Wide Model of Our NEXT CORE Series  
Available Now for Larger Users!




46cm Seat Width!  
Comfortable even  
for Larger Users!


Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled Front Seat Height 470 mm / 450 mm \* Seat cushion not included.  
NEXT-52B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥172,000 Caster x Rear Wheel [Inches] 6×24 (24×11/4)

Attendant-Propelled Front Seat Height 470 mm / 450 mm \* Seat cushion not included.  
NEXT-62B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥172,000 Caster x Rear Wheel [Inches] 6×16 (16×13/8)




Seat Upholstery Colors



Product Number	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width	Weight
NEXT-52B HB	460	400	445	250~310	480	365~455	960	1010	620	16.6kg
NEXT-62B HB	460	400	445	250~310	480	365~455	960	985	600	15.6kg

\* The factory setting of the front seat height is 470mm.  
\* Additional fees are required for adjustments to the front seat height before shipping.



Wonderful Features of NEXT CORE - Wide!

Position & Angle Adjustable Footrest Structure




The position and angle of the footrests are adjustable. Users can sit with their feet closer in towards themselves, or adjust the angle to find a position better for their joints.

Height Adjustable Armrest Structure



The armrest is adjustable according to the user's physical condition and the thickness of the cushion used. The armrest height is 4-level adjustable: 25 • 27 • 29 • 31 cm.

Puncture-Proof Tires with Excellent Shock Absorption



Our tires are more shock absorbant than conventional puncture-proof tires (high polymer tires). Hybrid tires are puncture-proof tires with the comfort level of air tires.

\* Refer to the chart of options available on P22.


NEXT-31B HB/31B (Self-Propelled)  
NEXT-41B HB/41B (Attendant-Propelled)

Flip-Up Armrests and Swing In-and-Away  
Footrests for Safer, Easier Transfers



The back frame largely curves towards the back to optimize seat comfort and ease of use.

NEXT-31B (Self-Propelled) Seat Upholstery Colors



Swing In-and-Away



Can be swung straight up and even detached!

Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled Front Seat Height 430 mm \* Seat cushion not included.  
NEXT-31B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥142,000 Caster x Rear Wheel [Inches] 6×22 (22×11/4)

Attendant-Propelled Front Seat Height 430 mm \* Seat cushion not included.  
NEXT-41B HB Manufacturer's Suggested Retail Price (Non-Taxable) ¥142,000 Caster x Rear Wheel [Inches] 6×16 (16×13/8)

Duralock Design (Foot Brakes)

Attendant-Propelled Front Seat Height 430 mm \* Seat cushion not included.  
NEXT-41B HBD Manufacturer's Suggested Retail Price (Non-Taxable) ¥162,000 Caster x Rear Wheel [Inches] 6×16 (16×13/8)

Air Tire Design

Self-Propelled Front Seat Height 430 mm \* Seat cushion not included.  
NEXT-31B Manufacturer's Suggested Retail Price (Non-Taxable) ¥136,000 Caster x Rear Wheel [Inches] 6×22 (22×11/4)

Attendant-Propelled Front Seat Height 430 mm \* Seat cushion not included.  
NEXT-41B Manufacturer's Suggested Retail Price (Non-Taxable) ¥136,000 Caster x Rear Wheel [Inches] 6×16 (16×13/8)

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-31B HB	430	380/400/420	380	400	235	440	320~410	880	950	560
NEXT-31B										550
NEXT-41B HB	430	380/400/420	380	400	235	440	320~410	880	945	560
NEXT-41B										550

\* Footrest length: distance between the footrests and seat  
© Specifications are subject to change without prior notice for continual product improvement.  
© The dimensions are based on the seat width set at 400mm.



\* Refer to the chart of options available on P22.



NEXT-41B  
(Attendant-Propelled)



NEXT-11B HB/11B (Self-Propelled)

NEXT-21B HB/21B (Attendant-Propelled)




\* May not be JIS compatible with optional parts installed.



## A Standard Model of NEXT CORE Series Both User-Friendly and Functional

### Seat Upholstery Colors



Polyester

### Newly Designed Wheels & Push Rings



The wheels are easier to maintain thanks to the fewer number of spokes. The newly designed hand rim is the result of painstaking effort for a much easier grip.

### Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled	Front Seat Height <b>430 mm</b>	* Seat cushion not included.	[Folded] 290	Weight 12.4kg
NEXT-11B HB	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥116,000</b>	Caster x Rear Wheel [Inches] <b>6 × 22 (22 × 11/4)</b>	655	

Attendant-Propelled	Front Seat Height <b>430 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 11.4kg
NEXT-21B HB	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥116,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	655	

### Duralock Design (Foot Brakes)

Attendant-Propelled	Front Seat Height <b>430 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 12.7kg
NEXT-21B HBD	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥136,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	655	

### Air Tire Design

Self-Propelled	Front Seat Height <b>430 mm</b>	* Seat cushion not included.	[Folded] 290	Weight 11.8kg
NEXT-11B	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥110,000</b>	Caster x Rear Wheel [Inches] <b>6 × 22 (22 × 11/4)</b>	655	

Attendant-Propelled	Front Seat Height <b>430 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 10.8kg
NEXT-21B	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥110,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	655	

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-11B HB	430	380/400/420	380	400	235	440	320~410	880	935	550
NEXT-11B										
NEXT-21B HB	430	380/400/420	380	400	235	440	320~410	880	935	540
NEXT-21B HBD										
NEXT-21B										

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

© The dimensions are based on the seat width set at 400mm.

\* Refer to the chart of options available on P22.

Curved Back Frame & Adjustable Backrest  
(Standard Feature in NEXT CORE Series)



The back frame largely curves towards the back to optimize seat comfort and ease of use.

### Adjustable Footrests

The wheelchair features newly designed footrests adjustable without the use of any tools.



### Slim & Lightweight

With a total width size of 550 mm, the wheelchair weight only 11.8kg\*.

We upgraded its parts and the shape of its frame to achieve a highly practical wheelchair.

\* For NEXT-11B (Self-Propelled) with a seat width of 400mm



**NEXT-21B**  
(Attendant-Propelled)

NEXT-12B (Self-Propelled) • NEXT-22B (Attendant-Propelled)



\* May not be JIS compatible with optional parts installed.



NEXT-12B (Self-Propelled)

## A High-Seat Standard Wheelchair with Quality Seating Comfort for Larger Users

### Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled	Front Seat Height <b>470 mm</b>	* Seat cushion not included.	[Folded] 300	Weight 13.4kg
NEXT-12B HB	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥137,000</b>	Caster x Rear Wheel [Inches] <b>6 × 24 (24 × 11/4)</b>	715	

Attendant-Propelled	Front Seat Height <b>470 mm</b>	* Seat cushion not included.	[Folded] 290	Weight 12.4kg
NEXT-22B HB	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥137,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	715	

### Air Tire Design

Self-Propelled	Front Seat Height <b>470 mm</b>	* Seat cushion not included.	[Folded] 300	Weight 12.8kg
NEXT-12B	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥131,000</b>	Caster x Rear Wheel [Inches] <b>6 × 24 (24 × 11/4)</b>	715	

Attendant-Propelled	Front Seat Height <b>470 mm</b>	* Seat cushion not included.	[Folded] 290	Weight 11.8kg
NEXT-22B	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥131,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	715	

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-12B HB	470	420	400	445	250	460	350~440	940	965	570
NEXT-12B										
NEXT-22B HB	470	420	400	445	250	460	350~440	940	940	560
NEXT-22B										

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

\* Refer to the chart of options available on P22.

NEXT-10B (Self-Propelled) • NEXT-20B (Attendant-Propelled)



\* May not be JIS compatible with optional parts installed.



NEXT-10B (Self-Propelled)

## A Low-Seat Standard Wheelchair with Shallow Seat Depth for Shorter Users

### Hybrid Tire Design (Puncture-Proof Tires)

Self-Propelled	Front Seat Height <b>390 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 12.4kg
NEXT-10B HB	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥137,000</b>	Caster x Rear Wheel [Inches] <b>6 × 20 (20 × 11/4)</b>	615	

Attendant-Propelled	Front Seat Height <b>390 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 11.7kg
NEXT-20B HB	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥137,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	615	

### Air Tire Design

Self-Propelled	Front Seat Height <b>390 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 11.8kg
NEXT-10B	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥131,000</b>	Caster x Rear Wheel [Inches] <b>6 × 20 (20 × 11/4)</b>	615	

Attendant-Propelled	Front Seat Height <b>390 mm</b>	* Seat cushion not included.	[Folded] 280	Weight 11.1kg
NEXT-20B	Manufacturer's Suggested Retail Price (Non-Taxable) <b>¥131,000</b>	Caster x Rear Wheel [Inches] <b>6 × 16 (16 × 13/8)</b>	615	

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-10B HB	390	400	360	365	235	440	325~415	860	925	550
NEXT-10B										
NEXT-20B HB	390	400	360	365	235	440	325~415	860	925	540
NEXT-20B										

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

\* Refer to the chart of options available on P22.



NEXT-71B (Self-Propelled)•NEXT-81B (Attendant-Propelled)  
NEXT-70B (Self-Propelled)•NEXT-80B (Attendant-Propelled)




Watch a Video  
of the Features

## An Incredibly Tight-Turning Wheelchair (Total Length 83cm\*)!! Spins in Place in Hallways Only 90cm Wide!!

\* When in Compact Mode

Seat Upholstery Colors



### Hybrid Tire Design (Puncture-Proof Tires)

[Mid-Seat]

Self-Propelled Front Seat Height **430 mm** \* Seat cushion not included.  
NEXT-71B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥170,000** 4 × 22 (22 × 11/4) × 4 [Folded] 340 Weight 19.6kg Caster x Rear Wheel [Inches] 675

Attendant-Propelled Front Seat Height **430 mm** \* Seat cushion not included.  
NEXT-81B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥170,000** 4 × 16 (16 × 13/8) × 4 [Folded] 340 Weight 18.6kg Caster x Rear Wheel [Inches] 675

[Low-Seat]

Self-Propelled Front Seat Height **405 mm** \* Seat cushion not included.  
NEXT-70B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥170,000** 4 × 20 (20 × 11/4) × 4 [Folded] 340 Weight 19.1kg Caster x Rear Wheel [Inches] 650

Attendant-Propelled Front Seat Height **400 mm** \* Seat cushion not included.  
NEXT-80B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥170,000** 4 × 12 (12 × 13/8) × 4 [Folded] 340 Weight 17.8kg Caster x Rear Wheel [Inches] 645

### Air Tire Design

[Mid-Seat]

Self-Propelled Front Seat Height **430 mm** \* Seat cushion not included.  
NEXT-71B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥164,000** 4 × 22 (22 × 11/4) × 4 [Folded] 340 Weight 19.0kg Caster x Rear Wheel [Inches] 675

Attendant-Propelled Front Seat Height **430 mm** \* Seat cushion not included.  
NEXT-81B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥164,000** 4 × 16 (16 × 13/8) × 4 [Folded] 340 Weight 18.0kg Caster x Rear Wheel [Inches] 675

[Low-Seat]

Self-Propelled Front Seat Height **405 mm** \* Seat cushion not included.  
NEXT-70B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥164,000** 4 × 20 (20 × 11/4) × 4 [Folded] 340 Weight 18.5kg Caster x Rear Wheel [Inches] 650

Attendant-Propelled Front Seat Height **400 mm** \* Seat cushion not included.  
NEXT-80B HB Manufacturer's Suggested Retail Price (Non-Taxable) **¥164,000** 4 × 12 (12 × 13/8) × 4 [Folded] 340 Weight 17.2kg Caster x Rear Wheel [Inches] 645

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
NEXT-71B HB	430	400	360	405	225~285	460	295~370	975	860	580
NEXT-71B										
NEXT-81B HB	430	400	360	405	225~285	460	295~370	975	860	580
NEXT-81B										
NEXT-70B HB	405	400	360	380	225~285	460	295~340	950	860	580
NEXT-70B										
NEXT-80B HB	400	400	360	375	225~285	460	295~340	945	860	580
NEXT-80B										

© Specifications are subject to change without prior notice for continual product improvement.  
© Alterations or custom orders are unavailable for this model.

Lead a More Active Life with NEXT CORE - Spin



The footrests are about 12cm closer towards the user than a conventional wheelchair to prevent the risk of buttocks sliding forward when putting the user's feet up onto the footrests. It is easier to keep the pelvis untitled, and to find a good position for eating or being active.

## Choosing NEXT CORE Based on Height

Users About 150cm Tall

Low-Seat Standard Wheelchair

NEXT CORE - Petit



NEXT-10B NEXT-20B

Manufacturer's Suggested Retail Price

¥131,000 (Non-Taxable)

Front Seat Height [mm]	Front Seat Height [mm]
390	390
Seat Depth [mm]	Seat Depth [mm]
360	360
Total Width [mm]	Total Width [mm]
550	540
Weight [kg]	Weight [kg]
11.8	11.1

Users About 160cm Tall

Standard Wheelchair

NEXT CORE



NEXT-11B NEXT-21B

Manufacturer's Suggested Retail Price

¥110,000 (Non-Taxable)

Front Seat Height [mm]	Front Seat Height [mm]
430	430
Seat Depth [mm]	Seat Depth [mm]
380	380
Total Width [mm]	Total Width [mm]
550	540
Weight [kg]	Weight [kg]
11.8	10.8

Users About 170cm Tall

High-Seat Standard Wheelchair

NEXT CORE - Grand




NEXT-12B NEXT-22B

Manufacturer's Suggested Retail Price

¥131,000 (Non-Taxable)

Front Seat Height [mm]	Front Seat Height [mm]
470	470
Seat Depth [mm]	Seat Depth [mm]
400	400
Total Width [mm]	Total Width [mm]
570	560
Weight [kg]	Weight [kg]
12.8	11.8

## Choosing NEXT CORE Based on Height (Rough Estimate) & Features






# 3 Tips for Sitting Stability on NEXT CORE

Independence Begins with a Stable Posture.




## A Frame That Doesn't Touch the Shoulders! Our Backward Curving Backrest Frame

The backrest frame curves largely backwards to prevent the frame from coming into contact with the user's shoulders. Combining this with backrest adjustment features allows even those with rounded backs to sit as far back into their seats as possible while keeping their backs straight.






## Supported from Both Sides! Pelvic & Core Stability

The gap between the user and the wheelchair is considerably reduced by the frame and seat. The added lateral support along with the backrest and seat stabilizes the posture of users who tend to lean towards one side.



## Independent Left-Right Backrest Adjustment! Supports Good Posture from 3 Directions with Our 3D Sling Seat


Straight on, from the side and the top. Each NEXT CORE wheelchair is equipped with a 3D Sling Seat backrest. The left and right sides of the backrest are independently adjustable. The backrest stabilizes the posture of users with rounded backs or scoliosis.



# How to Adjust the Backrest Upholstery (Recommendation)

\* The method explained below is just one example.  
We recommend adjusting the backrest with the cooperation of a medical care professional.

## 1. Set Up Before Adjusting the Backrest Upholstery



- Belt①: Adjust along the neutral line  
\* Set slightly loose for users with curved backs.
- Belt②: Fasten the belt on both sides of the seat.
- Belt③: Fasten the belt on both sides of the seat.
- Belt④: Adjust along the neutral line.



Tutorial videos on adjusting the backrest available on YouTube.


- 2. Sit straight back in the seat. The buttocks should be as far back in the seat as possible. Adjust the footrest length.  
(Unfastening the belts during the initial settings provides enough give for any back curvatures.)

- 3. Fasten belt ①.  
Decide on the approximate position of the user's head.

- 7. Readjust belts ③ and ④.  
Readjust belt ① if necessary.



- 6. Lightly support the lower part of the rib cage with both hands. Fasten belt ② diagonally downward to replace the support from each hand.



- 4. To prevent the pelvis from tilting backwards, fasten belt ③ in the direction of the pelvis (posterior superior iliac spine).



- 5. Fasten belt ④ diagonally downward to ensure the seat is a good fit for the buttocks.





## NEXT CORE Series Option Chart

	Hybrid Tires (Puncture-Proof)	Walking Cane Holder	Rear IV Pole Set	Oxygen Tank Holder	Anti-Tip Bars (Fall Prevention) with Casters	Brake Extension Bar (5cm • 10cm)	Mudguards (Made of Resin)	Spoke Cover (Clear Only)
NEXT-10B	●	●	●	●	●	●	●	×
NEXT-11B	●	●	●	●	●	●	●	●
NEXT-12B	●	●	●	●	●	●	●	×
NEXT-20B	●	●	●	●	●	●	●	×
NEXT-21B	●	●	●	●	●	●	●	×
NEXT-22B	●	●	●	●	●	●	●	×
NEXT-31B	●	●	●	●	●	●	×	●
NEXT-41B	●	●	●	●	●	●	×	×
NEXT-50B	●	●	●	●	●	●	×	×
NEXT-51B	●	●	●	●	●	●	×	●
NEXT-52B HB	Standard Issue	●	●	●	●	●	×	×
NEXT-60B	●	●	●	●	●	●	×	×
NEXT-61B	●	●	●	●	●	●	×	×
NEXT-62B HB	Standard Issue	●	●	●	●	●	×	×
NEXT-70B	●	×	×	×	×	×	×	×
NEXT-71B	●	×	×	×	×	×	×	●
NEXT-80B	●	×	×	×	×	×	×	×
NEXT-81B	●	×	×	×	×	×	×	×

© Specifications are subject to change without prior notice for continual product improvement.

\* We recommend installing anti-tip bars (fall prevention) with the oxygen tank holder. Please be careful of the wheelchair toppling backwards.

### Hybrid Tires (Puncture-Proof)



This option makes pumping air into the tires unnecessary. These puncture-proof airless tires are practically as shock absorbent as air tires. The inner tubes are highly durable with a service life of about 5 years (excluding tire wear).

\* The tire in the photo is different from that of NEXT CORE wheelchairs.

### NEXT CORE Walking Cane Holder (Large)



### Rear IV Pole Set



### For NEXT CORE Oxygen Tank Holder



\* We recommend installing anti-tip bars (fall prevention) with the oxygen tank holder. Please be careful of the wheelchair toppling backwards.

\* Foldable oxygen tank holders are also available.

### For NEXT CORE Anti-tip Bars (Fall Prevention) with Casters



### Brake Extension Bar (5cm • 10cm)



\* 10cm in the photo

### For NEXT CORE Spoke Cover (22-Inch, Clear Only)



### Seat Belts for 3D Sling Seats



Image of Installation on Chest (For Illustrative Purposes Only)

\* Please avoid using the seat belt for restraining purposes.

### For NEXT CORE Mudguards Self-Propelled • Attendant-Propelled (Made of Resin)



## Standard · Multi-Functional Wheelchairs





Steel Standard Model

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



DM-81

Front Seat Height 470 mm

DM-81

Manufacturer's Suggested Retail Price (Non-Taxable) ¥57,600

Caster x Rear Wheel [Inches] 7 × 24 (24 × 13/8)

[Folded] 330  
Weight 17.6kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
DM-81	470	420	400	440	235	390	365~	885	1,045	625

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.


Please note that alterations are unavailable for this model.



B-16

Feature of DM Series

Vinyl Leather Seat



The seat upholstery is made of vinyl leather. This is the part of the wheelchair that tends to get the dirtiest. It is easy to wipe the seat clean.

Options

Hybrid Tires (Puncture-Proof)



Equipped with special shock absorbent inner tubes, this option makes pumping air into the tires unnecessary.

\* Material: Polyurethane Material (PU Material)

For 1 Wheelchair

Suggested Retail Price (Non-Taxable)  
¥6,000 UP

Solid Tires



These puncture-proof tires are made of rubber. They are suitable for use in hospitals and care facilities.

\* The added weight is 1.5kg per wheelchair.

No additional fees for DM Series only.

Oxygen Tank Holder (Standard)



This metal fixture is for holding an oxygen tank on the rear side of a wheelchair. It is suitable for oxygen tanks with a radius of 110mm.

For 1 Wheelchair

Suggested Retail Price (Non-Taxable)  
¥10,800 UP

Anterior IV Pole Set



①Anterior IV Pole Set (1 Iron Pole)  
For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥7,200 UP

②Anterior IV Pole Set (Aluminum, Retractable)  
For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥11,200 UP

Posterior IV Pole Set



①Posterior IV Pole Set (1 Iron Pole)  
For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥9,000 UP

②Posterior IV Pole Set (Aluminum, Retractable)  
For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥13,000 UP

\* Aluminum retractable IV pole in the photo.

Anti-tip Bars (Fall Prevention)

Anti-tip Bars (Fall Prevention) with Rubber Tips



This metal anti-tip bars (fall prevention) fixture is retrofit onto wheelchairs to prevent them from toppling over. It is suitable for users who tend to shift their center of gravity backwards when going uphill.

For 1 Wheelchair Suggested Retail Price (Non-Taxable) ¥4,800 UP

Anti-tip Bars (Fall Prevention) with Casters



This anti-tip bars (fall prevention) pipe is equipped with casters. It is easy to move forward even when it is touching the ground.

For 1 Wheelchair Suggested Retail Price (Non-Taxable) ¥7,000 UP

Anti-tip Bars (Fall Prevention) Clamp Model



This model of anti-tip bars (fall prevention) is suitable for use on wheelchairs installed with oxygen tank holders and rear IV poles.


For 1 Wheelchair Suggested Retail Price (Non-Taxable) ¥10,000 UP

Custom Printed Backrest • Personalized Name Plate



Names can be printed onto the back of the backrest. It is easy to determine the owner of the wheelchair, contributing to theft prevention.

\* The photo is just one example.



A personalized name plate is attached to the sideguard. \* The photo is just one example.

Custom Printed Backrest  
For 1 Wheelchair Suggested Retail Price (Non-Taxable)  
From ¥3,000 UP

Personalized Name Plate  
For 1 Wheelchair Suggested Retail Price (Non-Taxable)  
From ¥3,000 UP



Flip-Up Armrests & Swing-Away Footrests


Self-Propelled



Front Seat Height **470 mm**

**CM-251HB** Manufacturer's Suggested Retail Price (Non-Taxable) **¥93,700**

Caster x Rear Wheel [Inches]  
**7 × 24 (24 × 13/8)**



Weight **19.9kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
CM-251HB	470	420	400	440	245	385	340~	880	1,025	625

\* Footrest length: distance between the footrests and seat  
© Specifications are subject to change without prior notice for continual product improvement.  
There is more weight in the back when installing a posterior IV pole or oxygen tank to the wheelchair.  
We recommend installing anti-tip bars (fall prevention).



B-36

Flip-Up Armrests & Swing-Away Footrests

Attendant-Propelled



Front Seat Height **470 mm**

**CM-261HB** Manufacturer's Suggested Retail Price (Non-Taxable) **¥93,700**

Caster x Rear Wheel [Inches]  
**7 × 16 (16 × 13/8)**



Weight **18kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
CM-261HB	470	420	400	440	245	385	340~	880	980	575

\* Footrest length: distance between the footrests and seat  
© Specifications are subject to change without prior notice for continual product improvement.  
There is more weight in the back when installing a posterior IV pole or oxygen tank to the wheelchair.  
We recommend installing anti-tip bars (fall prevention).



B-36

Flip-Up Armrests & Elevating Unit


Self-Propelled



Front Seat Height **470 mm**

**CM-251HB-SE** Manufacturer's Suggested Retail Price (Non-Taxable) **¥109,700**

Caster x Rear Wheel [Inches]  
**7 × 24 (24 × 13/8)**



Weight **21.7kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
CM-251HB-SE	470	420	400	440	245	385	360~	880	1,040	625

\* Footrest length: distance between the footrests and seat  
© Specifications are subject to change without prior notice for continual product improvement.  
There is more weight in the back when installing a posterior IV pole or oxygen tank to the wheelchair.  
We recommend installing anti-tip bars (fall prevention).



B-36

Flip-Up Armrests & Elevating Unit

Attendant-Propelled



Front Seat Height **470 mm**

**CM-261HB-SE** Manufacturer's Suggested Retail Price (Non-Taxable) **¥109,700**

Caster x Rear Wheel [Inches]  
**7 × 16 (16 × 13/8)**



Weight **19.8kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
CM-261HB-SE	470	420	400	440	245	385	360~	880	995	575

\* Footrest length: distance between the footrests and seat  
© Specifications are subject to change without prior notice for continual product improvement.  
There is more weight in the back when installing a posterior IV pole or oxygen tank to the wheelchair.  
We recommend installing anti-tip bars (fall prevention).



B-36

Features of Models CM-251HB and 261HB

Flip-Up Armrest



It is easier to transfer from the wheelchair to beds or toilets by flipping up the armrest. It alleviates strain on the attendant.


Hybrid Tires (Puncture-Proof)



Equipped with special shock absorbent inner tubes, this option makes pumping air into the tires unnecessary.

\* Material: Polyurethane Material (PU Material)

Swing-Away & Detachable Footrests



The footrests can be swung away to shift the wheelchair closer to beds and toilets for easier transfers. This feature is also useful for foot propelling.

(CM-251HB, 261HB)

Options for CM-251 and 261

Posterior IV Pole Set



①Posterior IV Pole Set (1 Iron Pole)  
②Posterior IV Pole Set (Aluminum, Retractable)

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)

**¥9,000 UP**

\* Aluminum retractable IV pole in the photo.

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)

**¥13,000 UP**

Flip-Up Oxygen Tank Holder



This metal fixture is for holding an oxygen tank on the rear side of a wheelchair. It is suitable for oxygen tanks with a radius of 110mm.

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable) **¥10,800 UP**

\* Standard oxygen tank holder in the photo.

Fall Prevention

Anti-tip Bars (Fall Prevention) with Rubber Tips

This metal fall prevention fixture is retrofit onto wheelchairs to prevent them from toppling over. It is suitable for users who tend to shift their center of gravity backwards when going uphill.

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable) **¥4,800 UP**



Anti-tip Bars (Fall Prevention) with Casters

This model of anti-tip bars is equipped with casters. It is easy to move forward even when touching the ground.

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable) **¥7,000 UP**



Anti-tip Bars (Fall Prevention) Clamp Model

This model of anti-tip bars is suitable for use on wheelchairs installed with oxygen tank holders and rear IV poles.

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable) **¥10,000 UP**





Aluminum Standard Model

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



AR-101

Front Seat Height 430 mm

AR-101

Manufacturer's Suggested Retail Price (Non-Taxable) ¥60,000

Caster x Rear Wheel [Inches] 6 x 22 (22 x 13/8)



Weight 13.8kg

Front Seat Height 380 mm

AR-111

Manufacturer's Suggested Retail Price (Non-Taxable) ¥60,000

Caster x Rear Wheel [Inches] 6 x 22 (22 x 13/8)



Weight 13.6kg

Unit: mm


Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-101	430	400/420	400	400	250	430	360~	890	1,015	645
AR-111	380	400/420	380	360	250	430	280~	855	960	645

\* Footrest length: distance between the footrests and seat.

\* The dimensions are based on the seat width set at 400mm.

© Specifications are subject to change without prior notice for continual product improvement.

Please note that alterations are unavailable for this model.



S-1

S-2

Options

Hybrid Tires (Puncture-Proof)



Equipped with special shock absorbent inner tubes, this option makes pumping air into the tires unnecessary.  
\* Material: Polyurethane Material (PU Material)

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥6,000 UP

Oxygen Tank Holder




This metal fixture is for holding an oxygen tank on the rear side of a wheelchair. It is suitable for oxygen tanks with a radius of 110mm.  
\* Foldable oxygen tank holders are also available.

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥10,800 UP

Custom Printed Backrest • Personalized Name Plate



Names can be printed onto the back of the backrest. It is easy to determine the owner of the wheelchair, contributing to theft prevention.  
\* The photo is just one example.



A personalized name plate is attached to the sideguard. \* The photo is just one example.

Custom Printed Backrest

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
From ¥3,000 UP

Personalized Name Plate

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
From ¥3,000 UP

Standard Wheelchair for Children



MKD-M

Aluminum Standard Model

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



AR-201B

Front Seat Height 430 mm

AR-201B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥78,000

Caster x Rear Wheel [Inches] 6 x 22 (22 x 13/8)



Weight 15.2kg

Front Seat Height 380 mm

AR-211B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥78,000

Caster x Rear Wheel [Inches] 6 x 22 (22 x 13/8)



Weight 14.7kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-201B	430	380/400/420	400	400	250	420	360~	900	1,015	670
AR-211B	380	380/400/420	380	360	250	425	285~	860	955	670

\* Footrest length: distance between the footrests and seat.

\* The dimensions are based on the seat width set at 400mm.

© Specifications are subject to change without prior notice for continual product improvement.

Please note that alterations are unavailable for this model.



S-1

S-2

High-Seat Model

Self-Propelled




Seat Upholstery Color / S-2

Front Seat Height 470 mm

AR-271B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥78,000

Caster x Rear Wheel [Inches] 7 x 24 (24 x 13/8)



Weight 16.0kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-271B	470	400/420/440	420	440	240	460	350~	930	1,020	690

\* The dimensions are based on the seat width set at 420mm.

\* Footrest length: distance between the footrests and seat.

\* The footrest length is at the shortest setting available.

© Specifications are subject to change without prior notice for continual product improvement.

Please note that alterations are unavailable for this model.



S-1

S-2

Options


Hybrid Tires (Puncture-Proof)



Equipped with special shock absorbent inner tubes, this option makes pumping air into the tires unnecessary.  
\* Material: Polyurethane Material (PU Material)

For 1 Wheelchair  
Suggested Price (Non-Taxable)  
¥6,000 UP

Walking Cane Holders



Walking Cane Holder (Large)

This is an excellent size for holding 2 crutches or 2 walking canes.



Walking Cane Holder (Small)

The is an excellent size for holding 1 walking cane.



Aluminum Standard Wheelchair AR Series

Aluminum Standard Model

Attendant-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



AR-301

Seat Upholstery Color / S-2




S-1 S-2

Front Seat Height 430 mm

AR-301

Manufacturer's Suggested Retail Price (Non-Taxable) ¥75,600

Caster x Rear Wheel [Inches] 6 × 16 (16 × 13/8)



Weight 13.8kg

Front Seat Height 380 mm

AR-311

Manufacturer's Suggested Retail Price (Non-Taxable) ¥75,600

Caster x Rear Wheel [Inches] 6 × 16 (16 × 13/8)



Weight 13.5kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-301	430	380/400/420	400	400	250	420	360~	895	1,015	600
AR-311	380	380/400/420	380	360	250	425	285~	860	940	600

Unit: mm

\* Footrest length: distance between the footrests and seat.  
\* The dimensions are based on the seat width set at 400mm.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.

High-Seat Model

Attendant-Propelled



Seat Upholstery Color / S-1



S-1 S-2

Front Seat Height 470 mm

AR-371

Manufacturer's Suggested Retail Price (Non-Taxable) ¥75,600

Caster x Rear Wheel [Inches] 7 × 16 (16 × 13/8)



Weight 14.8kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-371	470	400/420/440	420	440	240	460	350~	930	1,000	615

Unit: mm

\* The dimensions are based on the seat width set at 420mm.  
\* Footrest length: distance between the footrests and seat.  
\* The footrest length is at the shortest setting available.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.

Options

Hybrid Tires (Puncture-Proof)



Equipped with special shock absorbent inner tubes, this option makes pumping air into the tires unnecessary.

\* Material: Polyurethane Material (PU Material)

For 1 Wheelchair  
Suggested Retail Price (Non-Taxable)  
¥6,000 UP

Standard Feature for AR Series (Except AR-101-111)

Drumbrakes



This standard feature is for riding safely on slopes and streets wet from the rain.

Feature of AR-311

Shallow Seat Depth



The shallow seat depth is ideal for shorter users to sit as far back as possible in the seat and to maintain a more stable sitting posture.

(AR-311)

Aluminum Wide Wheelchair AR Series

Aluminum Wide Model

Self-Propelled



Load Capacity 130kg



Seat Upholstery Color / S-2



S-1 S-2

Front Seat Height 430 mm

AR-280

Manufacturer's Suggested Retail Price (Non-Taxable) ¥114,000

Caster x Rear Wheel [Inches] 6 × 22 (22 × 13/8)




Weight 16.3kg

Front Seat Height 380 mm

AR-301

Manufacturer's Suggested Retail Price (Non-Taxable) ¥75,600

Caster x Rear Wheel [Inches] 6 × 16 (16 × 13/8)



Weight 13.5kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-280	430	440/460/480	400	400	250	420	350~	885	1,010	725

Unit: mm

\* Footrest length: distance between the footrests and seat.  
\* The dimensions are based on the seat width set at 460mm.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.

Aluminum Wide Model

Attendant-Propelled



Seat Upholstery Color / S-1




S-1 S-2

Front Seat Height 430 mm

AR-380

Manufacturer's Suggested Retail Price (Non-Taxable) ¥111,600

Caster x Rear Wheel [Inches] 6 × 16 (16 × 13/8)



Weight 14.8kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-380	430	440/460/480	400	400	250	420	350~	885	1,010	660

Unit: mm

\* Footrest length: distance between the footrests and seat.  
\* The dimensions are based on the seat width set at 460mm.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.

Features of Models AR-501, 511B and 601

Flip-Up & Detachable Armrests



It is easier to transfer from the wheelchair to beds or toilets by flipping up the armrest towards the rear. It is also detachable to alleviate strain on the attendant.

(AR-501, 601)

Swing-Away & Detachable Footrests



The footrests can be swung away to shift the wheelchair closer to beds and toilets for easier transfers. This feature is also useful for foot propelling.

Shallow Seat Depth



The shallow seat depth is ideal for shorter users to sit as far back as possible in the seat and to maintain a more stable sitting posture.

(AR-511B)



Aluminum Multi-Functional Wheelchair AR Series

Flip-Up Armrests & Swing-Away Footrests

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



Seat Upholstery Color / S-1



Front Seat Height 420 mm

AR-501

Manufacturer's Suggested Retail Price (Non-Taxable) ¥110,400

Caster x Rear Wheel [Inches] 6 × 22 (22 × 13/8)



Weight 16.6kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-501	420	380/400/420	400	390	250	375	325~	850	955	625

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.  
When installing a rear IV pole or oxygen tank to the wheelchair, please use anti-tip bars (fall prevention).

Low-Seat Model

Self-Propelled



Seat Upholstery Color / S-1




Front Seat Height 380 mm

AR-511B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥92,400

Caster x Rear Wheel [Inches] 6 × 22 (22 × 13/8)



Weight 15.3kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-511B	380	380/400/420	380	360	235	430	260~	910	950	670

\* The dimensions are based on the seat width set at 400mm.  
\* Footrest length: distance between the footrests and seat  
\* The footrest length is at the shortest setting available.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.

Flip-Up Armrests & Swing-Away Footrests

Attendant-Propelled




\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



Seat Upholstery Color / S-2



Front Seat Height 420 mm

AR-601

Manufacturer's Suggested Retail Price (Non-Taxable) ¥105,600

Caster x Rear Wheel [Inches] 6 × 16 (16 × 13/8)



Weight 15.2kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-601	420	380/400/420	400	390	250	375	325~	850	945	550

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.  
When installing a rear IV pole or oxygen tank to the wheelchair, please use anti-tip bars (fall prevention).

Aluminum Semi-Modular Wheelchair AR Series

Semi-Modular Model

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 75kg



Seat Upholstery Color / S-1




Front Seat Height 420 mm

AR-901

Manufacturer's Suggested Retail Price (Non-Taxable) ¥117,600

Caster x Rear Wheel [Inches] 6 × 22 (22 × 13/8)




Weight 18.0kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-901	420	380/400/420	400	390	270~350	400	325~	850	955	610

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
\* The dimensions are subject to change depending the adjustment of each part.  
\* The armrest height is 5-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.  
Please note that alterations are unavailable for this model.  
When installing a rear IV pole or oxygen tank to the wheelchair, please use anti-tip bars (fall prevention).  
(For AR-901 Only)

3D Cushion (Standard Feature)




Size Adjustment Functions

Seat Width Adjustment

Adjustable by Switching to an X-Frame (Sold Separately)

The seat frame is fixed using e-Pins, so the seat width is adjustable (by switching to an X-frame) without using any tools. The tipping unit is specially designed into a form easy to step on.



Insert the pin inside the hole of the frame and rotate to lock securely.

Remove ❶. Pull ❷ out in the direction of the arrow.

Seat Height Adjustment

The adjustable caster fork is used to finely adjust the front seat height.



Set the caster installation spacer above or below the frame to roughly adjust the front seat height.



Backrest • Seat Upholstery Adjustment

The backrest and seat upholstery can be adjusted to fit each user's posture with separate belt buckles.



Height Adjustable Armrests

Armrests tend to be too low when using cushions with the wheelchair. This feature is useful for adjusting the armrests to a comfortable height. The armrests are made of high-strength resin for increased durability and precision.



Main Features of AR-901

Swing-Away & Detachable Footrests




The footrests can be conveniently swung away to shift the wheelchair closer to beds and toilets for much easier transfers. Remove the footrest on the healthy side of the user's body for a comfortable set-up for foot propelling.

Flip-Up Armrest



This features makes lateral transfers in and out of wheelchairs easier. The armrest is made of high-strength resin for increased durability and precision.

Detachable Armrest



Lateral transfers in and out of wheelchairs are made even easier for attendants by flipping up and detaching the armrest. (Mostly for two-person assist)

Drumbrakes



The drumbrakes are used by attendants to brake by controlling the wheel axle. This braking system is useful for riding safely on slopes and streets wet from the rain. Requiring little power input to brake, this feature is especially convenient for female and elderly attendants. (A standard feature on self- and attendant-propelled models)

Standard • Multi-Functional

Aluminum Semi-Modular Wheelchair




Aluminum Semi-Modular Wheelchair AR Series

Semi-Modular Model

Attendant-Propelled



Seat Upholstery Color / S-2




Front Seat Height **420 mm**

\* Seat cushion not included.

**AR-901**

Manufacturer's Suggested Retail Price (Non-Taxable) **¥117,600**

(Standard Settings)  
Caster x Rear Wheel [Inches] **6 × 16 (16 × 13/8)**




Weight **16.4kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-901	420	380/400/420	400	390	270~350	400	325~	850	945	550


\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
\* The dimensions are subject to change depending the adjustment of each part.  
\* The armrest height is 5-level adjustable in 20mm increments.  
\* Switching to an X-frame is necessary to adjust seat width.  
© Specifications are subject to change without prior notice for continual product improvement.  
**Please note that alterations are unavailable for this model.**  
When installing a rear IV pole or oxygen tank to the wheelchair, please use anti-tip bars (fall prevention). (For AR-901 Only)

3D Cushion (Standard Feature)




Main Features of AR-901

Installation Spacer Positions



Adjustable Fork



Axle Positions





Chart of Wheelchair Sizes (AR-901(Self-Propelled)(Attendant-Propelled))

	Front Seat Height	Rear Seat Height	Caster Size	Rear Wheel Size	Installation Spacer Position	Adjustable Fork Hole Position	Axle Position
Self-Propelled	390	365	5 Inches	20 Inches	A	a	1
	400	365	5 Inches	20 Inches	A	b	1
	410	385	6 Inches	20 Inches	A	b	2
	420	390	6 Inches	22 Inches	A	c	1
	430	410	6 Inches	22 Inches	A	d	2
	440	410	6 Inches	22 Inches	B	d	2
	450	420	6 Inches	24 Inches	B	e	1
	460	430	6 Inches	24 Inches	C	e	2
Attendant-Propelled	470	450	6 Inches	24 Inches	D	e	3
	390	360	5 Inches	16 Inches	A	a	3
	400	370	5 Inches	16 Inches	A	b	4
	410	380	6 Inches	16 Inches	A	b	4
	420	390	6 Inches	16 Inches	A	c	5
	430	400	6 Inches	16 Inches	A	d	5
	440	410	6 Inches	16 Inches	B	d	6
	450	420	6 Inches	16 Inches	B	e	6
	460	430	6 Inches	16 Inches	C	e	7
	470	440	6 Inches	16 Inches	D	e	7

\* The front seat height indicated refers to the height above the seat frame.

Options

Swing-Away & Elevating Unit A (Standard)

Swing-Away & Elevating Unit A (Standard) (Suggested Price ¥16,000 UP). This option can be adjusted to meet the needs of users with limited mobility in knees. Exchanging to swing-away footrests may require fine-tuning.



Flip-Up Oxygen Tank Holder



\* Standard oxygen tank holder in the photo.

Posterior IV Pole Set



\* Option/ Equipped with aluminum retractable IV pole

Anti-Tip Bars (Fall Prevention) (For AR-901 Only)



\* Those in the photo are equipped with casters.

Aluminum Compact Semi-Modular Wheelchair AR Series

Compact Semi-Modular Model

Self-Propelled



Seat Upholstery Color / S-2




Front Seat Height **370 mm**

\* Seat cushion not included.

**AR-911S**

Manufacturer's Suggested Retail Price (Non-Taxable) **¥117,600**

(Standard Settings)  
Caster x Rear Wheel [Inches] **5 × 18 (18 × 13/8)**



Weight **16.5kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-911S	370	380/400	360	350	200~280	455	270~	900	930	620

\* Footrest length: distance between the footrests and seat \* The dimensions are based on the seat width set at 400mm.  
\* The dimensions are subject to change depending the adjustment of each part.  
\* The armrest height is 5-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.  
**Please note that alterations are unavailable for this model.**  
When installing a rear IV pole or oxygen tank to the wheelchair, please use anti-tip bars (fall prevention). (For AR-911 Only)

Front Seat Height (mm)	Rear Seat Height (mm)	Caster Size	Installation Spacer Position	Adjustable Fork Position	Rear Wheel Size	Axle Position
330	310	5	A	*	18	①
350	330			C		②
370	350		C	E		③
390	370					④

3D Cushion (Standard Feature)



\* Will be equipped with a fork exclusively used for this size. Please refer to the chart below for information about A~B, a~e, ①~⑤.

Compact Semi-Modular Model

Attendant-Propelled



Seat Upholstery Color / S-1




Front Seat Height **370 mm**

\* Seat cushion not included.

**AR-911S**

Manufacturer's Suggested Retail Price (Non-Taxable) **¥117,600**

(Standard Settings)  
Caster x Rear Wheel [Inches] **5 × 16 (16 × 13/8)**



Weight **15.5kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-911S	370	380/400	360	350	200~280	455	270~	900	930	555

\* Footrest length: distance between the footrests and seat \* The dimensions are based on the seat width set at 400mm.  
\* The dimensions are subject to change depending the adjustment of each part.  
\* The armrest height is 5-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.  
**Please note that alterations are unavailable for this model.**  
When installing a rear IV pole or oxygen tank to the wheelchair, please use anti-tip bars (fall prevention). (For AR-911 Only)

Front Seat Height (mm)	Rear Seat Height (mm)	Caster Size	Installation Spacer Position	Adjustable Fork Position	Rear Wheel Size	Axle Position
330	310	5	A	*	16	②
350	330			C		③
370	350		E			④
390	370			⑤		

3D Cushion (Standard Feature)




\* Will be equipped with a fork exclusively used for this size. Please refer to the chart below for information about A~B, a~e, ①~⑤.

Seat Height Adjustment

The seat height of the wheelchair is adjustable by changing the size of the casters and rear wheels, or changing the mounting position. A self-propelled wheelchair can also be converted into a attendant-propelled one, and vice versa.


Assemblage of Installation Spacer



Axle Position




Adjustable Fork



Features of AR-911S

Extra-Low Seat

Shorter users can sit with their feet firmly planted in this extra-low seat wheelchair. It is safer and more comfortable to transfer in and out of.



Shallow Seat Depth

Shorter users can sit far back into the seat and maintain a more stable sitting posture.



Height Adjustable Armrests

The armrest height is adjustable to a low height (200mm) suitable for shorter users.



Standard • Multi-Functional

Aluminum Compact Semi-Modular Wheelchair



Stylish Multi-Functional Model

Self-Propelled

Seat Upholstery Color

M-2/M-3

M-1/M-3M-2/M-3

Front Seat Height **425 mm (445) mm**

\* Seat cushion not included.

**MY-1** Manufacturer's Suggested Retail Price (Non-Taxable) **¥140,800**

Caster x Rear Wheel [Inches] **6 (Super) × 22 (22 × 13/8)**

[Folded]  
(Front Seat Height 425)

370

Cushion Included

Weight **18.0kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MY-1	425 (445)	400	400	395/410	240~300	415	330~	835	1,065	650

\* Footrest length: distance between the footrests and seat.  
\* The armrest height is 20mm increments.  
\* The factory setting of the front seat height is 425mm.  
\* Additional fees are required for adjustments to the front seat height before shipping.  
© Specifications are subject to change without prior notice for continual product improvement.  
**There is more weight in the back when installing a posterior IV pole or oxygen tank to the wheelchair.**  
**We recommend installing anti-tip bars (fall prevention).**

3D Cushion (Standard Feature)

20 Rear

40 Front

Stylish Multi-Functional Model

Attendant-Propelled

Seat Upholstery Color

M-1/M-3

M-1/M-3M-2/M-3

Front Seat Height **425 mm (445) mm**

\* Seat cushion not included.

**MY-2** Manufacturer's Suggested Retail Price (Non-Taxable) **¥140,800**

Caster x Rear Wheel [Inches] **6 (Super) × 16 (16 × 13/8)**

[Folded]  
(Front Seat Height 425)

280

Cushion Included

Weight **16.3kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MY-2	425 (445)	400	400	395/410	240~300	415	330~	835	1,065	565

\* Footrest length: distance between the footrests and seat.  
\* The armrest height is 20mm increments.  
\* The factory setting of the front seat height is 425mm.  
\* Additional fees are required for adjustments to the front seat height before shipping.  
© Specifications are subject to change without prior notice for continual product improvement.  
**There is more weight in the back when installing a posterior IV pole or oxygen tank to the wheelchair.**  
**We recommend installing anti-tip bars (fall prevention).**

3D Cushion (Standard Feature)

20 Rear

40 Front

Features of MY Series

Adjustable Seat Height

The front seat height can be adjusted to 425/445mm.  
\* A standard issue wheelchair is set at 425mm.

Flip-Up Armrest

It is easier to transfer onto beds and toilets by flipping up the armrest. It alleviates strain on the attendant.

Swing-Away & Detachable Footrests

The footrests can be swung away to shift the wheelchair closer to beds and toilets for easier transfers. This feature is also useful for foot propelling.

Height Adjustable Armrests

The armrest can be adjusted to a height suitable for the user.

36

37



Lightweight Standard Wheelchair

Aerial Series



Lightweight Standard Model

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



Front Seat Height 430 mm

MW-SL11B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥127,200

Caster x Rear Wheel [Inches] 665 6×22 (22×1)

[Folded] 360  
Weight 11.0kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MW-SL11B	430	400	400	400	235	400	310~	840	975	650

\* Footrest length: distance between the footrests and seat  
\* The armrest height is 4-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.



E-2

Lightweight Standard Model

Attendant-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



Front Seat Height 430 mm

MW-SL21B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥121,200

Caster x Rear Wheel [Inches] 665 6×16 (16×1.75)

[Folded] 320  
Weight 10.0kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MW-SL21B	430	400	400	400	235	400	310~	840	975	595

\* Footrest length: distance between the footrests and seat  
\* The armrest height is 4-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.



E-2

Option

Option (For All Models of the Aerial Series)

● Adjustable Backrest Upholstery  
(Back Cushion · 3D Seat Cushion Included)

Weight .....0.5kg UP

Suggested Retail Price.....¥20,000 UP



Lightweight Multi-Functional Wheelchair

Aerial Series



Lightweight Flip-Up Armrest & Swing-Away Footrests

Self-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



Front Seat Height 430 mm

MW-SL31B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥163,200

Caster x Rear Wheel [Inches] 665 6×22 (22×1)

[Folded] 365  
Weight 12.2kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MW-SL31B	430	400	400	405	235	400	305~	840	975	580

\* Footrest length: distance between the footrests and seat  
\* The armrest height is 4-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.



E-2

Lightweight Flip-Up Armrest & Swing-Away Footrests

Attendant-Propelled



\* May not be JIS compatible with optional parts installed.

Load Capacity 100kg



Front Seat Height 430 mm

MW-SL41B

Manufacturer's Suggested Retail Price (Non-Taxable) ¥157,200

Caster x Rear Wheel [Inches] 665 6×16 (16×1.75)

[Folded] 310  
Weight 11.1kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MW-SL41B	430	400	400	405	235	400	305~	840	975	525


\* Footrest length: distance between the footrests and seat  
\* The armrest height is 4-level adjustable in 20mm increments.  
© Specifications are subject to change without prior notice for continual product improvement.



E-2

Features of Models MW-SL31B, SL41B, SL5B and SL6B

Flip-Up Armrest



It is easier to transfer from the wheelchair to beds or toilets by flipping up the armrest towards the rear. It is also detachable to alleviate strain on the attendant.

Swing-Away & Detachable Footrests



The footrests can be swung away to shift the wheelchair closer to beds and toilets for easier transfers. This feature is also useful for foot propelling.

Drumbrakes



This standard feature is for riding safely on slopes and streets wet from the rain.





## Lightweight & Extra-Low Seat Wheelchair with Flip-Up Armrest & Swing-Away Footrests

### Self-Propelled



Front Seat Height **370** mm **(350)** mm

**MW-SL5B** Manufacturer's Suggested Retail Price (Non-Taxable) **¥175,200**

Caster x Rear Wheel [Inches]  
**5 × 20 (20 × 1)**

[Folded]  
(Front Seat Height 370)

410

665

We  
12

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
<b>MW-SLSB</b>	370(350)	400	380	345/325	235~275	450	295~	840	930	610

- \* Footrest length: distance between the footrests and seat
- \* The armrest height is 5-level adjustable in 10mm increments.
- \* The footrest length is 4-level adjustable at 15mm increments.
- \* The factory setting of the front seat height is 370mm.
- \* Additional fees are required for adjustments to the front seat height before shipping.

© Specifications are subject to change without prior notice for continual product improvement.



E-2

## Lightweight & Extra-Low Seat Wheelchair with Flip-Up Armrest & Swing-Away Footrests

## Attendant Propelled



Front Seat Height **370** mm **(350)** mm

**MW-SL6B** Manufacturer's Suggested Retail Price (Non-Taxable) **¥169,200**

Caster x Rear Wheel [Inches]  
**5 × 16 (16 × 1.75)**

[Folded]  
(Front Seat Height 37")  
340  
665  
We  
11.

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
<b>MW-SL6B</b>	370(350)	400	380	345/325	235~275	450	295~	840	930	540

- \* Footrest length: distance between the footrests and seat
- \* The armrest height is 5-level adjustable in 10mm increments.
- \* The footrest length is 4-level adjustable at 15mm increments.
- \* The factory setting of the front seat height is 370mm.
- \* Additional fees are required for adjustments to the front seat height before shipping.

© Specifications are subject to change without prior notice for continual product improvement.



E-2

### Features of Models **MW-SL5B** and **SL6B**

### Height Adjustable Armrests

The armrest can be adjusted to a height suitable for the user. The armrest height can also be adjusted according to the thickness of the cushion used.  
(235mm~275mm/5 levels of adjustment)



## Invertable Pads


The direction of armrest pads can be inverted from left to right, and vice versa. The armrest can be used as a support for users to stand up independently by extending it towards the frontside of the wheelchair.




## Choosing a Wheelchair Made Super Easy!

# Matsunaga's Guide to Choosing a Wheelchair


## A—Standard Model Self-Propelled NEXT CORE Series



**B—Standard Model** Attendant-Propelled NEXT CORE Series



## C—Tilt & Reclining Model OASIS Series & My Tilt 3D Series





# TOPICS



GOOD DESIGN AWARD 2016  
**GOLD AWARD**

(PIIF) Onogi Science and Technology Promotion Foundation  
Winner of the Top Silver Award in 2016



**GERMAN  
DESIGN  
AWARD  
WINNER  
2018**



## Bringing Japanese Hospitality to the World.

From airport counters to  
onboard aircrafts,  
our barrier-free wheelchair  
morph  
made of 100% resin is  
the way to go.

Product information  
available on YouTube.






## Reclining & Tilt Wheelchair





My Tilt - Mini 3D



MH-SRL, MH-SRL-SE Only

Equipped with Lift-Up Unit

This feature enables lateral movements. The wheelchair can be placed right next to the bed and rotate in place.

Easy to Use Even for Shorter Users  
Suitable User Height **140cm~165cm**

\* Rough Estimate

Front Seat Height **410 mm** \* Seat cushion not included  
**MH-SR** Manufacturer's Suggested Retail Price (Non-Taxable) **¥198,000** Caster x Rear Wheel [Inches] **5 × 12 (12 × 13/8)** [Folded] 340 990

Front Seat Height **410 mm** \* Seat cushion not included  
**MH-SR-SE** Manufacturer's Suggested Retail Price (Non-Taxable) **¥214,000** Caster x Rear Wheel [Inches] **5 × 12 (12 × 13/8)** [Folded] 340 990

Front Seat Height **410 mm** \* Seat cushion not included  
**MH-SRL** Manufacturer's Suggested Retail Price (Non-Taxable) **¥218,000** Caster x Rear Wheel [Inches] **5 × 12 (12 × 13/8)** [Folded] 340 990

Front Seat Height **410 mm** \* Seat cushion not included  
**MH-SRL-SE** Manufacturer's Suggested Retail Price (Non-Taxable) **¥234,000** Caster x Rear Wheel [Inches] **5 × 12 (12 × 13/8)** [Folded] 340 990

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width	Weight
MH-SR	410	400	350	410	0-140~300	480	260~	1,010	890	530	24.7kg
MH-SR-SE	410	400	350	410	0-140~300	480	300~	1,010	880~1,195	530	26.5kg
MH-SRL	410	400	350	410	0-140~300	480	260~	1,010	890	530	27.0kg
MH-SRL-SE	410	400	350	410	0-140~300	480	300~	1,010	880~1,195	530	29.0kg


Reclining Angle	Tilt Angle
90°~110°	0°~25°



\* The Swing-Away & Elevating Unit B (90 Degree) for exclusive use with My Tilt - Mini 3D is also sold individually.  
\* Please note that alterations are unavailable for this model.

Easy to Use Even for Shorter Users! Adjustable Features of My Tilt - Mini 3D!

● A Seat Shorter Users Can Sit Far Back In Shallow Seat Depth (35cm)



[Deeper Seat] The back of the knees touch the seat. The user cannot sit far back.  
[Shallower Seat] The user can sit far back to reduce the risk of sacral sitting (slumping).

● Tucked in as Much as Possible 90 Degree Footrests (Position • Angle Adjustable)



Users unable to extend their knees or suffering from malformation of feet joints (equinus contracture) can properly position their feet on the footrests.  
\* The footrests can be slid 5cm forward from the position closest to the user.



Conventional Product Minus 10.5cm

● Freely Adjustable Back & Forth • Left & Right from a Low Position Snake Headrest

The headrest is adjustable to a low position optimal for shorter users. Thanks to this new structure, the position and angle headrest can be adjusted left & right • back & forth simultaneously.




● Supports Stable Sitting Positions for Users with Rounded Back Forward-Folding Backrest Frame (Positive Feature)



The Forward-Folding Backrest Frame is a popular mechanism in our OASIS - Positive Series. It is easier to reduce the gap between the backrest and the shoulder girdle which tend to form for users with rounded backs. This stabilizes the head to neck areas and makes it easier to find an optimal eating position.

● Hybrid Tire Design  
\* Puncture-Proof Tires (Standard Feature)

My Tilt - Compact 3D



Front Seat Height **430 mm (450) mm** \* Seat cushion not included.  
**MH-CR3D** Manufacturer's Suggested Retail Price (Non-Taxable) **¥192,000** Caster x Rear Wheel [Inches] **6 × 12 (12 × 13/8)** [Folded] 325 835 Weight **24.5kg**

\* The factory setting of front and rear seat height is 430mm.

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MH-CR3D	430mm (450)	400	380	430-450	0,140~300	560	330~	1,085	1,020	530

\* Footrest length: distance between the footrests and seat  
\* The armrest height is 8-level adjustable in 20mm increments.  
\* The initial setting of the total length is 1,020mm.  
\* The factory setting of the front seat height is 430mm.  
\* Additional fees are required for adjustments to the front seat height before shipping.  
© Specifications are subject to change without prior notice for continual product improvement.

Reclining Angle	Tilt Angle
90°~125°	0°~30°

3D Cushion (Standard Feature)



1. Left & Right Side of Backrest Independently Adjustable
2. Pelvis & Core Support Frame
3. Backward Curving Backrest Frame

Duralock • Mini Tackle Brake (Standard Feature)

It is possible to brake the wheelchair from both the front and back.

Swing-Away Footrests (Standard Feature)

Up until now, swing-away footrests were an option on My Tilt wheelchairs. Now they are a standard feature.

\* The options of swing-away & elevating units can also be installed.

Tilt & Reclining Feature



30° Tilt



30° Tilt  
130° Recline

Equipped with features the equivalent of current My Tilt wheelchair for similar sitting stability.

Fold for an Even More Compact Wheelchair!



The wheelchair can be loaded onto kei cars (small cars). It makes going out e.g. to the hospital more convenient.

Options



■ Swing-Away Footrests & Elevating Unit A (Standard)  
This unit can be adjusted for users with limited movement in their knees.  
Suggested Retail Price (Non-Taxable) **¥16,000up**  
\* The added weight is about 2.0kg.  
\* If this option is selected at the time of purchase, the standard swing-away feature is not included.



■ Hybrid Tire (Puncture-Proof)  
This option makes pumping air into the tires unnecessary. These puncture-proof airless tires are practically as shock absorbent as air tires.  
Suggested Retail Price (Non-Taxable) **¥6,000up**  
\* The added weight is about 0.4kg.



■ Oxygen Tank Holder for MH-CR3D  
\* Also Compatible with MH-SR Series  
Suggested Retail Price **Discretionary Price**



■ EXPOLE  
Suggested Retail Price **Discretionary Price**



■ Snake Headrest  
Suggested Retail Price (Non-Taxable) **¥40,000up**



■ Seat Belts for 3D Sling Seats  
Suggested Retail Price **Discretionary Price**  
\* Please avoid using the seat belt for restraining purposes.




# My Tilt - Grand 3D

## My Tilt - Grand 3D




A Tilt & Reclining Wheelchair Taller Users Can Sit Comfortably In  
Equipped with a lift-up unit for lateral movement and in-place rotation.



Front Seat Height 470 mm	* Seat cushion not included	Swing-Away Design
<b>MH-GR</b>	Manufacturer's Suggested Retail Price (Non-Taxable)	¥199,000
Caster x Rear Wheel [Inches]		6×16 (16×13/8)
Front Seat Height 470 mm	* Seat cushion not included	Swing-Away Design & Elevating Unit
<b>MH-GR-SE</b>	Manufacturer's Suggested Retail Price (Non-Taxable)	¥215,000
Caster x Rear Wheel [Inches]		6×16 (16×13/8)
Front Seat Height 470 mm	* Seat cushion not included	Swing-Away Design & Lift-Up Unit
<b>MH-GRL</b>	Manufacturer's Suggested Retail Price (Non-Taxable)	¥235,000
Caster x Rear Wheel [Inches]		6×16 (16×13/8)
Front Seat Height 470 mm	* Seat cushion not included	Swing-Away & Elevating Design • Lift-Up Unit
<b>MH-GRL-SE</b>	Manufacturer's Suggested Retail Price (Non-Taxable)	¥251,000
Caster x Rear Wheel [Inches]		6×16 (16×13/8)

Product Number	Seat Height	Seat Width	Seat Depth	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width	Tilt Angle	Reclining Angle	Weight
MH-GR	470	440	380*	0-140~300	630	380~	1115	1030	620	0°~30°	95°~130°	25.0kg
MH-GR-SE	470	440	380*	0-140~300	630	330~	1115	980	620	0°~30°	95°~130°	26.7kg
MH-GRL	470	440	380*	0-140~300	630	380~	1115	1030	620	0°~30°	95°~130°	29.5kg
MH-GRL-SE	470	440	380*	0-140~300	630	330~	1115	980	620	0°~30°	95°~130°	31.2kg

\* The seat cushion is 470mm deep.



### Super User-Friendly Features on My Tilt - Grand!

#### Lift-Up Unit

\* MH-GRL, MH-GRL-SE Only

**Lift-Up Bar**

When the lift-up bar is lowered, the wheels float up.

This enables lateral movement and in-place rotation!

#### Back & Forth, Left & Right Integrated Footbrake

**Pinch**

The brake can be activated · deactivated by stepping on any of the back · front · left · right brakes.

#### One-Touch Height Adjustable BJ Headrest

**Pinch**

One touch is all it takes to adjust the headrest height. It is very easy to adjust the headrest position after transferring onto the wheelchair.

#### Slope Rollers

Even when the underside of the footrests bump into the slope, the slope rollers rotate to help to keep moving up the slope. They lower the risk of damages.

### Options

**Snake Headrest**

Suggested Price (Non-Taxable) **¥40,000 UP**

**Oxygen Tank Holder for MH-GR**

\* Incompatible with MH-GRL, MH-GRL-SE Series

Suggested Price **Discretionary Price**

## 3D Model 3 Tips for Sitting Stability

### 1 Backward Curving Backrest Frame

The backrest frame curves largely backwards to prevent the frame from coming into contact with the shoulders.

Conventional Backrest Frame      My Tilt - Compact 3D

The user can sit comfortably by reducing the pressure on the shoulder blades. This alleviates pressure to the abdomen.

### 2 Core Support Frame

This support wraps itself around the user's pelvis to hips from both sides. It stabilizes the sitting position of users who tend to incline to one side.

Image of Core Support Frame (For Illustrative Purposes Only)

### 3 3D Support for Proper Sitting Positions! 3D Sling Seat

\* The 3D Sling Seat is a registered trademark of Deai-Kobo, LLC.


**【Comparing the Pressure from Conventional Seats and 3D Sling Seats】**  
The backrest properly supports the entire back to reduce pressure on the buttocks.

Belts can be fastened only where the user needs support

Belts can be loosened only where they need to be

Conventional Seat      3D Sling Seat


## My Tilt Series: Sizing Guide





My Tilt - Rigid 3D

My Tilt - Rigid 3D



Front Seat Height **450 mm (430-470) mm**

\* Seat cushion not included

**MH-RD3D** Manufacturer's Suggested Retail Price (Non-Taxable) **¥189,000** Caster x Rear Wheel (Inches) **4 × 16 (16 × 13/8)**

Product Number	Front Seat Height	Seat Width	Seat Depth	Armrest Height	Backrest Height	Foot Length	Total Height	Total Length	Total Width	Weight
MH-RD3D	450 (430-470)	420	380 (350)	0-140~300	560	280~	1,105 (Front Seat Height 450)	940	600	24.1kg

Unit: mm

\* The factory setting of the front seat height is 450mm.  
\* Additional fees are required for adjustments to the front seat height before shipping.  
© Specifications are subject to change without prior notice for continual product improvement.  
© Alterations or custom orders are unavailable for this model.


■ Tilt Angle 0°~30° ■ Reclining Angle 90°~130°

①~④ Features of My Tilt - Rigid!


Adjustable Seat Height  
(Factory Setting at 450mm)



Adjustable Seat Depth  
(Factory Setting at 380mm)



One-Touch Height Adjustable BJ Headrest



One touch is all it takes to adjust the headrest height. It is very easy to adjust the position of the headrest after transferring onto the wheelchair.

Back & Forth, Left & Right Interlocking Footbrakes




The brake can be activated · deactivated by stepping on any of the back · front · left · right brakes.

90 Degree Footrests



\* Swing-Away Design Only  
Even users unable to extend their legs can position their feet on these footrests.

3D Model



- Conventional Backrest Frame      My Tilt - Compact 3D
- 1 Left & Right Sides of the Backrest Independently Adjustable
  - 2 Pelvis & Core Support Frame
  - 3 Backward Curving Backrest Frame

Plenty of User-Friendly Options



■ Snake Headrest  
Suggested Retail Price (Non-Taxable)  
**¥40,000 up**



■ Oxygen Tank Holder For MH-RD3D  
Suggested Retail Price (Non-Taxable)  
**Discretionary Price**



■ Swing-Away & Elevating Unit A (Standard)  
Suggested Retail Price (Non-Taxable)  
**¥16,000 up**



■ EXPOLE  
Suggested Retail Price (Non-Taxable)  
**Discretionary Price**



■ Seat Belts for 3D Sling Seats  
Suggested Retail Price (Non-Taxable)  
**Discretionary Price**  
\* Please avoid using the seat belt for restraining purposes.

\* Please contact us to learn more about other options.

Tilt & Reclining Wheelchair My Tilt - Value

My Tilt - Value



Front Seat Height **430 mm**

\* Seat cushion not included

**MH-VR** Manufacturer's Suggested Retail Price (Non-Taxable) **¥178,000** Caster x Rear Wheel (Inches) **6 × 16 (16 × 13/8)**

Front Seat Height **430 mm**

\* Seat cushion not included

**MH-VR-SE** Manufacturer's Suggested Retail Price (Non-Taxable) **¥194,000** Caster x Rear Wheel (Inches) **6 × 16 (16 × 13/8)**

Unit: mm


Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width	Weight
MH-VR	430	420	400	425	0-180~240	475	370~	1,150	1,140	605	23.6kg
MH-VR-SE	430	420	400	425	0-180~240	475	330~	1,150	1,135	605	25.0kg

■ Tilt Angle 0°~30° ■ Reclining Angle 90°~125° ■ Seat Cushion Thickness 30mm

© Please note that alterations are unavailable for this model.



The angle is adjusted based on the user's physical condition.  
Sliding forward can be prevented.  
\* 30° tilt, 125° recline in the photo.




The elevating feature allows the footrests to be adjusted for users with limited movement in their knees.  
\* A standard feature of MH-VR-SE.

Valuable Features on a Budget. Our Original My Tilt Wheelchair with Highly-Demanded Standard Features.



● Hybrid Tire (Puncture-Proof)  
Pumping air into these puncture-proof tires is unnecessary. They are always safe to ride on.



● Back & Forth Interlocking Footbrakes  
The brake can be easily activated · deactivated just by stepping on any of the back · front · left · right brakes without crouching down.



● BJ Headrests  
Freely adjust the angle and its left · right · back & forth · up & down position of the headrest according to the position of the user's head.  
\* Pillow Shape: BJ- B Type (Small Concave) Standard



■ Oxygen tank holder for MH-VR  
Suggested Price  
**Discretionary Price**

● Drumbrakes  
This brake is for stable control even in sloped areas.

● Adjustable Backrest Upholstery  
The backrest upholstery can be adjusted by fastening belts according to the shape of the user's back.  
\* We recommend users with curved and uneven backs our wheelchairs with 3D sling seats for independently adjustable left and right sides.



## 1 Interlocking Tilt & Reclining System

Reclining only engages the pelvic area to tilt to prevent the user from sliding forward. Minimal friction occurs between the backrest and back when reclining. The user will more likely maintain good posture even after repeated reclining with this mechanism.



## 2 Choice of Headrest Shapes

Choose between two BJ-A and BJ-B headrests.

BJ-A  
(Large Concave)

The headrest is largely concave. The head can be positioned far back for good stability.



BJ-B  
(Small Concave)

The headrest features a shallow concave. It is easier for the user to move the head.




## 3 Height Adjustable Headrest Pipe

Choose between our straight and curved-back models.

Straight Model



Curved-Back Model  
(80mm Bend)



## 4 BJ Headrest

The angle and back & forth position of the headrest can be adjusted. Position the headrest for a sitting position with the user's chin tucked in to reduce the risk of aspiration.



## 5 Swing-Away & Detachable Footrests

The footrests can be swung away to shift the wheelchair closer to beds and toilets for easier transfers. This feature is also useful for foot propelling.




## 6 Adjustable Backrest Upholstery

The backrest upholstery can be adjusted to fit the user's posture for more comfortable seating.

## 7 Angle & Width Adjustable Footrests

Slide the footrests left and right to ensure no gaps form between each footrest. This is especially useful for modular wheelchairs with adjustable seat widths.




## 8 Forward-Folding Joint on the Upper Section of the Backrest (OS-11TRSP · 12TRSP)

Equipped with our newly designed forward-folding joint, the upper section of the backrest is curved forward to enable the entire backrest to support an active forward-leaning posture. Combining this with OASIS double reclining has remarkably improved its body pressure distribution performance.



## 9 Equipped with Footbrakes (Duralock) [Attendant-Propelled Model Only]

The drumbrake is coupled with the footbrake for more stable braking control.



## 10 Ultra-Compact

The total width of the Attendant-Propelled model is 530mm. That of the self-propelled model is 570mm. Thanks to the compact design, both models are excellent for home use.



## Compact Wheelchair with Interlocking Tilt & Reclining System

OASIS Series

### Forward-Folding Joint Model

Self-Propelled



Seat Upholstery Color

V-1 Headrest / BJ-B Height Adjustable Pipe / Straight

\* Please choose between two types of headrests (BJ-A or BJ-B) and two types of height adjustable pipes (straight and curved back).



Front Seat Height <b>390 mm (370) mm</b>		Caster x Rear Wheel [Inches]		915		Weight	
OS-11TRSP		Manufacturer's Suggested Retail Price (Non-Taxable)		¥234,600		23.5kg	


Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
OS-11TRSP	390mm (370)	380/400/420	350	390	180~240	550	285~	1,025	925	615

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
\* The armrest height is 4-level adjustable at 20mm increments.  
\* The headrest is not included in the total height of the wheelchair.  
\* The seat cushion height is 50mm.  
\* The factory setting of the front seat height is 390mm.  
\* Additional fees are required for adjustments to the front seat height before shipping.  
© Specifications are subject to change without prior notice for continual product improvement.

Interlocking	
Reclining Angle	Tilt Angle
90°~115°	0°~25°

### Height Adjustable Headrest Pipe

Straight Model



Curved-Back Model (80mm Bend)



### Shape of Headrest

BJ-A (Large Concave)



BJ-B (Small Concave)



• Waveform Plastic Push Ring

### Seat Upholstery Color

Option +¥3,000		
N-2 (Navy) Knit	V-1 (Red) Polyester	B-61 (Navy) Vinyl

### Forward-Folding Joint

Attendant-Propelled



Seat Upholstery Color

N-2 Headrest / BJ-B Height Adjustable Pipe / Straight

\* Please choose between two types of headrests (BJ-A or BJ-B) and two types of height adjustable pipes (straight or curved back).



Front Seat Height <b>430 mm</b>		Caster x Rear Wheel [Inches]		955		Weight	
OS-12TRSP		Manufacturer's Suggested Retail Price (Non-Taxable)		¥212,600		22.7kg	


Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
OS-12TRSP	430	380/400/420	350	430	0-140~260	550	300~	1,065	905	530

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
\* The armrest height is 7-level adjustable at 20mm increments.  
\* The headrest is not included in the total height of the wheelchair.  
\* The seat cushion height is 50mm.  
© Specifications are subject to change without prior notice for continual product improvement.


Interlocking	
Reclining Angle	Tilt Angle
90°~115°	0°~25°

### Height Adjustable Headrest Pipe

Straight Model



Curved-Back Model (80mm Bend)



### Shape of Headrest

BJ-A (Large Concave)



BJ-B (Small Concave)



### Shallow Seat Depth

Shorter users can sit far back into the seat for a more stable sitting position.

### Seat Upholstery Color

Option +¥3,000		
N-2 (Navy) Knit	V-1 (Red) Polyester	B-61 (Navy) Vinyl

\* Refer to the chart of options available on p.53.



Compact Wheelchair with Interlocking Tilt & Reclining System

OASIS Series



Compact Model with Interlocking Tilt & Reclining System


Attendant-Propelled



Seat Upholstery Color/N:2 Headrest/BJ-A Height Adjustable Pipe/Straight

\* Please choose between two types of headrests (BJ-A or BJ-B) and two types of height adjustable pipes (straight or curved back).

Height Adjustable Headrest Pipe



Shape of Headrest

BJ-A (Large Concave) BJ-B (Small Concave)



Seat Upholstery Color



Front Seat Height **430 mm** \* Seat cushion not included  
**OS-12TR** Manufacturer's Suggested Retail Price (Non-Taxable) **¥182,600** Caster x Rear Wheel [Inches] **6 × 12 (12 × 13/8)** [Folded] 320 Weight 22.6kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
OS-12TR	430	380/400/420	400	430	0-140~	450	330~	965	995	530

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
\* The armrest height is 6-level adjustable at 20mm increments.  
\* The headrest is not included in the total height of the wheelchair.  
\* The seat cushion height is 50mm.  
© Specifications are subject to change without prior notice for continual product improvement.

Interlocking	
Reclining Angle	Tilt Angle
90°~115°	0°~25°

Compact Size of Total Width 530mm




This ultra-slim reclining wheelchair has a total width of 530mm. It is easy to operate in narrow spaces. Remove the seat cushion and headrest when folding the wheelchair up for storage.

\* Refer to the chart of options available on p.53.

Standard Features of the OASIS Series

New Mechanism to Prevent Forward Sliding & Relieve Abdominal Pressure



Reclining only engages the pelvic area to tilt to prevent the user from sliding forward. Minimal friction occurs between the backrest and back when reclining. The user will more likely maintain good posture even after repeated reclining with this mechanism.

Adjust the Headrest Angle & Position for Better Eating Posture



The angle and back & forth position of the headrest can be adjusted. Position the headrest for a sitting position with the user's chin tucked in to reduce the risk of aspiration.

Options

OASIS Series



OASIS Reclining Options (Same for All Models in Series)

Option: Seat Upholstery



●B-61 Vinyl Suggested Retail Price (Non-Taxable) **¥3,000 UP**



●V-1 Polyester Suggested Retail Price (Non-Taxable) **¥3,000 UP**

Snake Headrest



Suggested Retail Price (Non-Taxable) **¥40,000 UP**

Swing-Away & Elevating Unit



Suggested Retail Price (Non-Taxable) **¥16,000 UP**

\* The added weight is about 2.0kg.  
\* If this option is selected at the time of purchase, the standard swing-away feature is not included.

Mini Tackle Brake

[Attendant-Propelled Wheelchairs Only]



Suggested Retail Price (Non-Taxable) **¥10,000 UP**

Hybrid Tires (Puncture-Proof Tires)



Suggested Retail Price (Non-Taxable) **¥6,000 UP**

Table Set for OASIS Series



Suggested Retail Price (Non-Taxable) **¥18,800 UP**

Oxygen Tank Holder for OASIS Series

[Attendant-Propelled Wheelchairs Only]



Suggested Retail Price (Non-Taxable) **¥16,000 UP**

EXPOLE



Suggested Retail Price (Non-Taxable)  
**Discretionary Price**

\* Please make an inquiry to the distributor for more information about other options.  
\* The prices shown refer to the additional fees when purchasing a new wheelchair.



Steel Reclining Wheelchair CM Series

Reclining Model

Self-Propelled

Front Seat Height **470 mm**

Backrest · Footrests Move Separately

**CM-501HB** Manufacturer's Suggested Retail Price (Non-Taxable) **¥127,200**

Caster x Rear Wheel [Inches] **7 × 24 (24 × 13/8)**

[Folded]  
335

Weight  
26.4kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
CM-501HB	470	400	400	450	230	840	415~	1,270	1,155	655

Unit: mm

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

Reclining Angle

95°~160°

B-36

Reclining Model

Attendant-Propelled

Front Seat Height **470 mm**

Backrest · Footrests Move Separately

**CM-541HB** Manufacturer's Suggested Retail Price (Non-Taxable) **¥127,200**

Caster x Rear Wheel [Inches] **7 × 16 (16 × 13/8)**

[Folded]  
275

Weight  
24.5kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
CM-541HB	470	400	400	450	230	835	415~	1,270	1,090	585

Unit: mm

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

Reclining Angle

95°~160°

B-36

Aluminum Reclining Wheelchair MW Series

Reclining Model

Attendant-Propelled

Front Seat Height **460 mm**

Backrest · Footrests Move Separately

**MW-14HB** Manufacturer's Suggested Retail Price (Non-Taxable) **¥130,800**

Caster x Rear Wheel [Inches] **7 × 16 (16 × 13/8)**

[Folded]  
250

Weight  
18.5kg

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MW-14HB	460	400	400	435	240	850	405~	1,270	1,065	560

Unit: mm

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

Reclining Angle

90°~160°

Seat Upholstery Color / J-1

J-1 J-3

For All Models of the CM • MW Series

Hybrid Tires (Puncture-Proof Tires)

This option makes pumping air into the tires unnecessary. These puncture-proof airless tires are practically as shock absorbent as air tires. (Exchanges due to tire wear are not taken into account.)

Options (For All Models of the CM • MW Series)

Oxygen Tank Holder for Reclining Models

This oxygen tank holder can be retrofit onto reclining wheelchairs. Install onto the piping of the backrest frame using seat fixture screws.  
\* We recommend using anti-tip bars (fall prevention) when installing a oxygen tank holder.

Suggested Retail Price (Non-Taxable)

**¥15,600 UP**

(Anti-Tip Bars (Fall Prevention))

Suggested Retail Price (Non-Taxable)

**¥17,800 UP**

(Anti-Tip Bars (Fall Prevention) with Casters)

IV Pole Set for Elevating Models

This IV pole set (IV pole + bracket) can be retrofit onto elevating wheelchairs.  
\* Only anterior IV poles are available.  
\* Other devices can be installed for fall prevention.  
Please make an inquiry if interested.

Anterior IV Pole Bracket for Elevating Models

Suggested Retail Price (Non-Taxable)

**¥9,000 UP**

Anti-Tip Bars (Fall Prevention)

Anti-Tip Bars (Fall Prevention) with Rubber Tips

Suggested Retail Price (Non-Taxable)

**¥4,800 UP**

Anti-Tip Bars (Fall Prevention) with Casters

Suggested Retail Price (Non-Taxable)

**¥7,000 UP**

Reclining & Tilt Aluminum Reclining Wheelchair

54

55



Full Reclining Wheelchair Escort



Full Reclining Model

Attendant-Propelled



Seat Upholstery Color/B-36



B-36 Vinyl Leather

Front Seat Height **475 mm**

\* Seat cushion not included

**FR-11R**

Manufacturer's Suggested Retail Price (Non-Taxable) **¥162,000**

Caster x Rear Wheel [Inches] **6 × 12** (Puncture-Proof Tires)

Weight **29.0kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
FR-11R	475	400	420	455	0-180~240	825	335~	1,245	1,205	535

\* Footrest length: distance between the footrests and seat

\* The dimensions are based on the reclining angle at 100° on a cushion.

\* The armrest height is 5-level adjustable at 20mm increments.

\* The footrest length is 3-level adjustable at 20mm increments.

© Specifications are subject to change without prior notice for continual product improvement.



Reclining Angle

100°~170°

\* The total length is 1,800mm when the reclining angle is at 170°.

Features of Escort

Reclining Backrest Interlocked with Elevating Legrests



The backrest can be reclined smoothly without levels. The reclining angle can be adjusted to an almost fully flat mode. This model is recommended for heavily handicapped users on bedrest.

Hybrid Tires (Puncture-Proof Tires)



These puncture-proof tires are a standard feature. This option makes air pumping unnecessary.

5-Level Height Adjustable Armrests



The armrests can be lowered to the same height as the seat. Thanks to this feature, transferring in and out of the wheelchair is much easier. This alleviates strain on the attendants. (The armrest height is 5-level adjustable: 0mm, 180mm, 200mm, 220mm, and 240mm.)

Tilt & Reclining Wheelchair Capable of Lateral Movement Elise

Elise



Tilt & Full Reclining Model

Attendant-Propelled



Seat Upholstery Color/Elise Skin



Elise Skin Vinyl Leather (B-16)

Front Seat Height **475 mm**

\* Seat cushion not included

**FR-31TR**

Manufacturer's Suggested Retail Price (Non-Taxable) **¥237,600**

Caster x Rear Wheel [Inches] **6 × 12 (12 × 13/8)**

Weight **36.5kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
FR-31TR	475	400	400~420	465	0-40~240	850	335~	1,240	1,180	650

\* Footrest length: distance between the footrests and seat

\* The dimensions are based on the reclining angle at 100°, the tilt angle at 0°, on a cushion.

\* The armrest height is 10-level adjustable at 25mm increments.

\* The footrest length is 3-level adjustable at 20mm increments.

\* The factory setting of the seat depth is 410mm.

\* The total length when fully flat is 1,800mm.

© Specifications are subject to change without prior notice for continual product improvement.



Seat Upholstery Color/Vinyl Leather

Reclining Angle

100°~175°

Tilt Angle

0°~30°

Features of Elise

Tight turns are possible even in narrow spaces with this lift-up feature.

Normal



● Lower the lift-up bar.

Lift-Up



● The safety wheels lower to lift the rear wheels up.

This feature enables lateral movements. The wheelchair can be easily placed right next to the bed.



There is no need to repeatedly reposition the wheelchair when transferring onto the bed.

Tilt & Reclining Mechanism to Maintain a Relaxing Posture

● 45° Recline · 15° Tilt



● Full Recline



Footrests

The footrests can be folded backwards to free up space for assistance, such as when loading/unloading the wheelchair.



Armrests

The armrests can be lowered to the same height as the seat. Assistance is made easier. The armrest is 10-level adjustable to a height appropriate for the user's body type and posture. (0 · 40~240/25mm increments)



\* Range of reclining angle 100°~175°, range of tilt angle 0°~30°

Reclining & Tilt

Tilt & Reclining Wheelchair Capable of Lateral Movement



## NOTE

---







LUUK-01 (Self-Propelled) • LUUK-02 (Attendant-Propelled)

Matsunaga's Kids Series  
Modular Wheelchairs for Children

LUUK.  
is a great  
first  
wheelchair!

A modular  
frame  
that adjusts to  
growing bodies!

Just  
one frame  
to meet the  
needs of users  
90cm-130cm  
tall!

Wide Frame  
with Wide Footrests  
&  
Front & Back Adjustable  
Footrests

Support the Entire Back

Self-Propelled  
**LUUK-01** Manufacturer's Suggested  
Retail Price (Non-Taxable) **¥184,000~**

Attendant-Propelled  
**LUUK-02** Manufacturer's Suggested  
Retail Price (Non-Taxable) **¥184,000~**

Specification Chart Unit: mm

	Wheel Size	Front Seat Height	Armrest Height	Seat Width	Seat Depth	Backrest Height
Self-Propelled Model	4×20Inches	375 395	210	240	240	350
	4×22Inches	415 435	220	260	260	
	5×20Inches	455 475	230	280	280	415
	5×22Inches	410 420		300	300	
		440 460		320	320	
Attendant-Propelled	4×16Inches	400	240	340	320	480

\* Custom purchase orders are available.

Be Unique with Seat Upholstery Colors

● Standard Colors ● Optional Colors ¥10,000up

1. Black/Gray 2. Light Blue /Gray 3. Pink/Gray 4. Blue/Gray 5. Lime Green /Gray 6. Red/Gray

Be Unique with Frame Colors

● Standard Colors ● Optional Colors 1 ¥25,000up ● Optional Colors 2 ¥30,000up

1. Ivory 2. Shadow Black 3. Aqua Milk 4. Lime Green 5. Pearl Yellow 6. Peach Pink 7. Lavender

8. MP White 9. MP Gun Metallic 10. Sparkle Sky Blue 11. Sparkle Lime 12. Sparkle Orange 13. Candy Red

\* We have an assortment of soft pastel colors and vivid metallic colors.

Plenty of Options for New Wheelchairs ~List of Additional Fees~

 Drumbrake for Attendants Suggested Retail Price <b>¥16,000up</b>	 Footbrakes Suggested Retail Price <b>¥10,000up</b> *Attendant-propelled models come with duralock footbrakes. Note 1) Walking cane holders • type A anti-tip bars (fall prevention) cannot be installed	 Cushion Caster Suggested Retail Price <b>¥15,000up</b>	 Mudguards (Attendant-Propelled Model Only) Suggested Retail Price <b>¥10,000up</b>	 Type A Fall Prevention Suggested Retail Price <b>¥15,000up</b> Note) Walking cane holders and footbrakes cannot be installed
 Metal Fittings Set for Table Suggested Retail Price <b>¥14,000up</b> * Table top not included	 Chest Belt Suggested Retail Price <b>¥6,000up</b>	 Pelvic Belt Suggested Retail Price <b>¥6,000up</b>	 Foot Extension Prevention Belt Suggested Retail Price <b>¥7,000up</b>	 Seat Depth Lengthening Belt (40mm) Suggested Retail Price <b>¥10,000up</b>
Walking Cane Holder (Self-Propelled Only) ¥4,200 UP		Push Ring Coating ¥14,000 UP		High-Mount Footrests ¥10,000 UP
YAMAHA Spoke Cover				

\* Custom purchase orders of Luuk are available. Please inquire for more information.

1.2.3 Ways LUUK is User-Friendly!



Feature 1 Wide Frame Structure with Wide Footrests

Thanks to the wide frame structure with wide footrests, there is enough space for users to place both feet even with braces or knee separators. It is easier and safer to transfer into the wheelchair.

Feature 2 Safe & Secure Design with Less Risk of Falling.

The footrests are positioned closer to the user than the casters. Even when the user steps onto the footrests, the wheelchair does not easily topple over with this highly safe design.



Feature 3 3D Sling Seats and 90 Degree Footrests for Good Posture.

● 3D Sling Seat  
The shape of the frame and seat easily fit the curvature of the user's back to properly support the entire back.

● 90 Degree Footrests  
The footrests are positioned to allow users to sit comfortably with their feet tucked closer in.

Modular Features for Long-Term Use!

 Height Adjustable Front Caster Wheels	 Height • Front & Back Adjustable Rear Wheels

\* The above figure is a self-propelled wheelchair with 4-inch casters.



Lightweight Stroller for Children  
Stroller for Children with Posture Support Features

MB Series

MB Series PAT Registered



Lightweight Stroller



Seat Upholstery Color/N-1  
Frame Color/Pearl White



6.4kg  
Ultralight!!

Super Light! Folds with One Touch!  
A Very Convenient Stroller for Outings!

Front Seat Height 500 mm

MB-PONY

Manufacturer's Suggested Retail Price (Non-Taxable) ¥96,800

Caster x Rear Wheel [Inches] 5×5




Weight 6.4kg

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width	Backrest Angle
MB-PONY	500	300	230/260/290	400	560	190~300	900	950	420	90/95/100°

\* Footrest length: distance between the footrests and seat  
\* The factory setting of the seat depth is 230mm, and the backrest angle is 95°.  
\* The dimensions shown are measured from above the seat frame (pipe) and not from above the cushion.  
\* The thickness of the seat cushion is 10mm.  
© Specifications are subject to change without prior notice for continual product improvement.



Standard Features

- 3-Point Belt Included
- Adjustable Backrest Angle (90° · 95° · 100°)

Options

Sunshade



Suggested Retail Price (Non-Taxable)  
¥18,000UP

Bag



Suggested Retail Price (Non-Taxable)  
¥7,200UP

Pillow



Suggested Retail Price (Non-Taxable)  
¥9,000UP

Stroller for Children with Posture Support Features



Seat Upholstery Color/M-6

Front Seat Height 450 mm

NEW MB-1

Manufacturer's Suggested Retail Price (Non-Taxable) ¥90,200

Caster x Rear Wheel [Inches] 5×5





Weight 7.5kg

Unit: mm

Product Number	Front Seat Height	Sea Width	Seat Depth	Rear Seat Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width	Seat Angle
NEW MB-1	450	300	290/310/330	350	615	90~240	995	920	470	90/95/100°

\* Footrest length: distance between the footrests and seat  
\* The factory setting of the seat depth is 290mm, and the backrest angle is 95°.  
\* The dimensions shown are measured from above the seat frame (pipe) and not from above the cushion.  
\* The thickness of the seat cushion is 30mm.  
© Specifications are subject to change without prior notice for continual product improvement.



Reclining Model NEW MB-1R

Manufacturer's Suggested Retail Price (Non-Taxable) ¥121,000



Fix this reclining kit onto the armrest of a standard stroller.

\* The 4-point seat belt · super casters are options.



Handlebar

Use the handlebar located on the back section to control the stroller.



Reclining Angle

The backrest can be reclined smoothly without levels at 95°~122°.

\* Please note the kit can only be fitted at our company factory.

Tilting Stroller

Mignon

Tilting Stroller Mignon



Mignon  
Seat Upholstery Color/Rose  
Optional Accessories

Frame Color


Black only

Lightweight · Compact · Stylish Tilting Stroller with Great Posture Support Features.

Mignon S

Manufacturer's Suggested Retail Price (Non-Taxable) ¥149,000

Caster x Rear Wheel [Inches] 5×6




Weight 9.4kg

Mignon M

Manufacturer's Suggested Retail Price (Non-Taxable) ¥149,000

Caster x Rear Wheel [Inches] 5×6




Weight 10.0kg

Mignon L

Manufacturer's Suggested Retail Price (Non-Taxable) ¥149,000

Caster x Rear Wheel [Inches] 5×6




Weight 10.2kg

Main Features

- 1 The body support unit fits the user seamlessly.
- 2 The new seat structure gently wraps around the user to provide comfortable seating.
- 3 The wheelchair can be compactly folded for convenient transport and storage.

Detailed Diagram of Mignon

Load Capacity 40kg



New Comfortable Seat Structure

The seat structure consists of 3 layers: sweat-absorbing and fast-drying material, highly breathable cushion, and 3D mesh. Enjoy comfortable seating at all time with its high breathability and cushioning properties.



Standard Features



Chest · Pelvic Belt



Core · Pelvic Pad



Lumbar Parts



Anchor Parts

Options



Sunshade



Headrest



Shoulder Belt



Aspirator Mounting Stand (Foldable)  
Dimension:320×220×50mm, Load Capacity:30kg

Color Lineup



Size Chart

Size	Total Width	Total Length	Total Height	Seat Width	Seat Depth	Backrest Height	Footrest Height	Weight	Suitable User Height	Load Capacity
S	500mm	725mm	1,080mm	330mm	200mm~275mm	600mm	135mm~265mm	9.4kg	85cm~110cm	40kg
M	530mm	765mm	1,150mm	360mm	260mm~335mm	670mm	180mm~315mm	10.0kg	95cm~120cm	40kg
L	550mm	765mm	1,180mm	380mm	300mm~375mm	700mm	195mm~330mm	10.2kg	110cm~140cm	40kg

[Adjustable Parts]Adjustable backrest upholstery, position of core · pelvic pad, seat depth, footrest height, reclining angle (0°, 10°), tilt angle (~45° backward incline), legrest angle (85°, 90°), lumbar · anchor adjustment  
[Folded Dimensions]S:500×286×750mm M:530×288×778mm L:550×288×816mm

Child

Tilting Stroller

63

Child

Lightweight Stroller for Children / Stroller for Children with Posture Support Features

62



# Standard Wheelchair for Children

## Indoor 6-Wheel Drive Wheelchair for Children

Leveret  
Tight Turn Jr.

### Standard Model for Children

Self-Propelled




Seat Upholstery Color/AP-3  
Frame Color/Pastel Blue

Front Seat Height **400 mm**

**MKD-M**

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥152,400**

Caster x Rear Wheel [Inches] **5×20 (20×13/8)**



Weight  
**12.8kg**

Unit: mm


Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
<b>MKD-M</b>	400	300	305	380	200	350	200~250	900	815	515

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.



AP-3



Pastel Blue  
(MATS-073)



AP-4



Pastel Pink  
(MATS-074)

### Standard Model for Children

Self-Propelled




Seat Upholstery Color/AP-4  
Frame Color/Pastel Pink

Front Seat Height **425 mm**

**MKD-L**

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥152,400**

Caster x Rear Wheel [Inches] **6×22 (22×13/8)**



Weight  
**13.4kg**

Unit: mm


Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
<b>MKD-L</b>	425	320	330	405	200	350	250~300	925	835	535

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.



AP-3



Pastel Blue  
(MATS-073)



AP-4



Pastel Pink  
(MATS-074)

### Tight Turn Jr.



6-Wheel Drive for Tight Turns Indoors!




Seat Upholstery Color/AP-4  
Frame Color/Pastel Pink

Front Seat Height **380 mm**

**TT-Jr**

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥140,800**

Caster x Rear Wheel [Inches] **4×20 (20×1)**



Weight  
**12.0kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
<b>TT-Jr</b>	380	260/280/300	260/280/300	360	195	400	200~250	930	785	430

\* Footrest length: distance between the footrests and seat

\* The dimensions are based on the seat width at 260mm and the gap between the wheel and the push rings at 15mm.

\* The seat width and seat depth are interlinked. (E.g.: The seat depth is 260mm when the seat width is 260mm.)


© Specifications are subject to change without prior notice for continual product improvement.



AP-3



AP-4




Pastel Blue  
(MATS-073)



Pastel Pink  
(MATS-074)

### Main Features of Tight Turn Jr.

#### J Tackle Brake (For TT-Jr Only)



The brake shoes are placed further towards the back to prevent brake-shoe-related injuries while self propelling.

#### The wheelchair can be compactly folded.




Side View



Front View

\* The spoke cover is an option.

#### Minimum Turn Radius 510mm










Wheelchair with Automatic Braking Device

Auto Break Series

The Auto Break Series are an assistive device which prevents falls by automatically locking the rear wheels of the wheelchair even when the user abruptly stands up without engaging the brake lock. Please make sure to use the parking brake on the wheelchair when transferring on and out of the wheelchair.

Feature 1


The lock deactivates when sitting down!



The brake activates when standing up!


Feature 2

The wheelchair can be moved around without a user sitting on it. The automatic braking device can be deactivated using a lever.



Feature 3

Users can select the plate size.



\* S size will be installed if no preferences are indicated.

Aluminum Standard Model

Self-Propelled



S-1



S-2



Seat Upholstery Color/S-1

The cushion cover is changed to a vinyl fabric.

Example

Front Seat Height **430 mm**

**AR-201BT**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥114,000**

Caster x Rear Wheel [Inches] **6×22(22×13/8)**



Weight **18.7kg**

(The cushion is not included in the dimensions.) \* Cushion Thickness with Sensor 50mm

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-201BT	430	380/400/420	400	400	250	420	360~	900	1,015	670

\* Footrest length: distance between the footrests and seat

\* The dimensions are based on the seat width set at 400mm.

© Specifications are subject to change without prior notice for continual product improvement.

Aluminum Low-Seat Swing-Away Wheelchair

Self-Propelled



S-1



S-2



Seat Upholstery Color/S-1

The cushion cover is changed to a vinyl fabric.

Example


Front Seat Height **380 mm**

**AR-511BT**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥128,400**

Caster x Rear Wheel [Inches] **6×22(22×13/8)**



Weight **18.8kg**

(The cushion is not included in the dimensions.) \* Cushion Thickness with Sensor 50mm

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-511BT	380	380/400/420	380	360	210	430	260~	910	950	680

\* The dimensions are based on the seat width set at 400mm.

\* Footrest length: distance between the footrests and seat

\* The footrest length is set at the shortest possible setting.

\* The weight is that of a wheelchair fitted with plate size S.

© Specifications are subject to change without prior notice for continual product improvement.

Highly Durable & Safe Aluminum Wheelchair



S-1



S-2



Example

Front Seat Height **430 mm**

**AR-101BT**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥106,000**

Caster x Rear Wheel [Inches] **6×22(22×13/8)**



Weight **18.6kg**

(The cushion is not included in the dimensions.) \* Cushion Thickness with Sensor 50mm

Example


Front Seat Height **380 mm**

**AR-111BT**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥106,000**

Caster x Rear Wheel [Inches] **6×22(22×13/8)**



Weight **18.4kg**

(The cushion is not included in the dimensions.) \* Cushion Thickness with Sensor 50mm

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
AR-101BT	430	400/420	400	400	250	430	360~	890	1,015	675
AR-111BT	380	400/420	380	360	255	430	280~	855	960	675

\* Footrest length: distance between the footrests and seat

\* The dimensions are based on the seat width set at 400mm.

© Specifications are subject to change without prior notice for continual product improvement.

Extra Low-Seat Wheelchair for Shorter Users to Keep Feet Firmly Planted



E-2



Example

Front Seat Height **370 mm (350) mm**

**MW-SL5BT**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥211,200**

Caster x Rear Wheel [Inches] **5×20(20×1)**



Weight **16.1kg**

(The cushion is not included in the dimensions.) \* Cushion Thickness with Sensor 50mm

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MW-SL5BT	370(350)	400	380	345/325	235~	450	295~	840	930	625

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

\* The above wheelchair is only one example. This device can be fitted onto other wheelchair models. Please don't hesitate to inquire for more information.

3T (Tesla) MRI-Compatible Wheelchair

MRI Series

3T (Tesla) MRI-Compatible Wheelchair

Aluminum



This wheelchair is compatible with 3T (Tesla) MRI.



Seat Upholstery Color/B-6

Front Seat Height **430 mm**

**MR-3T**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥145,000**

Caster x Rear Wheel [Inches] **6×22(22×11/4)**



Weight **13.4kg**

Features of the 3T (Tesla) MRI Series

- Non-magnetic aluminum material is used in the main frame. It is safe to use in 3T (Tesla) MRI rooms.

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MR-3T	430	400	400	400	235	430	360~	890	1,015	650

\* Footrest length: distance between the footrests and seat

\* The footrest length is set at the shortest possible setting.

© Specifications are subject to change without prior notice for continual product improvement.

Please note that alterations are unavailable for this model.

3T (Tesla) MRI-Compatible Ultra High-Spec Wheelchair

Aluminum



This is an ultra high-spec model of 3T (Tesla) MRI wheelchairs.



Seat Upholstery Color/B-38

Front Seat Height **430 mm**

**MR-3TMG**

Manufacturer's Suggested Retail Price (Non-Taxable)

**¥395,000**

Caster x Rear Wheel [Inches] **6×22(22×1.75)**



Weight **13.9kg**

Features of the Ultra High-Spec 3T (Tesla) MRI Series

- The frame, bolt, nuts and other small parts are made of non-magnetic material. This minimizes the risk of MRI magnetic accidents even in 3T (Tesla) MRI rooms.

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
MR-3TMG	430	420	400	400	235	435	350	890	1,010	655

\* Footrest length: distance between the footrests and seat

© Specifications are subject to change without prior notice for continual product improvement.

Please note that alterations are unavailable for this model.

IV Pole for Ultra High-Spec 3T (Tesla) MRI Wheelchair



**MR-IG-3TMG**

(Tax Included) **¥165,000**

Weight **5.1kg**

- Frame Material/Aluminum
- Total Width 430 x Total Length 430 x Total Height 1,430~2,370mm
- Percentage of Non-Magnetic Material/100%



\* Stretchers are also available.

Special Features

3T (Tesla) MRI-Compatible Wheelchair

67


66




# One-Arm Drive Wheelchair

SCULL Type1

Aluminum



Lever Drive Product Approved Under Act for  
the Comprehensive Support of Persons with Disabilities




Front Seat Height **470 mm**

**0H-1** Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥198,000**


Caster x Rear Wheel [Inches]  
**6 × 22 (22 × 1)**

[Folded]  
325



Weight  
**16.3kg**

Seat Upholstery Color / M-4  
Frame Color / Pastel Orange



M-4

## Features of One-Arm Drive Scull

- Shift the lever in the direction the user wishes to advance to move back & forth or left & right.
- The wheelchair can be operated for short periods of training.
- To engage the brake while operating the wheelchair, shift the level in the direction opposite of the direction the user is advancing. The parking brake on both rear wheels can also be engaged using one hand.

Unit: mm

Product Number	Front Seat Height	Seat Width	Seat Depth	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Height	Total Length	Total Width
0H-1	470	380/400/420	400	450	240	435	325~	965	1,030	620

\* Footrest length: distance between the footrests and seat  
\* The dimensions are based on the seat width set at 400mm.  
\* The seat cushion thickness is not included in the dimensions. The seat cushion is 40mm thick.  
© Specifications are subject to change without prior notice for continual product improvement.

Special Features

One-Arm Drive Wheelchair

68







Foldable Model For Exterior Use



The Ultimate Walker for Preventative Care

**MV-100** Manufacturer's Suggested Retail Price (Non-Taxable) **¥42,000**


7 Inches

Weight 8.8kg

8 Inches

Weight 9.1kg

- Handgrip Height Adjustment**  
The handgrip height can be adjusted at 11-level increments from 750~890mm (8-inch model) to closely fit the user's body type and posture. Use the user's leg joint area as a reference point to adjust the height to ensure the walker is easy to push.
- Seat Height Adjustment**  
The seat can be adjusted to a height easy for the user to sit down and stand up. Adjust the seat height to 450mm or 470mm for the 8-inch model, or 440mm or 460mm for the 7-inch model.
- Special Bag Included**  
A bag with an original design is included as a standard feature to store small items. It is especially useful during small outings and shopping.



Product Number	Width Between Handles	Total Width	Total Length	Total Height	Seat Height
<b>8-Inch Model</b>	420	600	690	750~890 (11 Levels)	450/470 (2 Levels)
<b>7-Inch Model</b>	420	600	690	740~880 (11 Levels)	440/460 (2 Levels)

\* The total height changes according to the handgrip height (11-level adjustable at 15mm increments). \* Weight: The weight of the special bag is not included. (The bag weight 200g.)  
\* The factory setting of the seat height is 450mm (8-inch model) or 440mm (7-inch model). © Specifications are subject to change without prior notice for continual product improvement.

Compact Storage Under Eaves and in Front Entrances

The walker can be compactly folded to a total width of 270mm. It can be stored in narrow spaces like front entrances and under eaves.



Seat Upholstery Color/M-4

Easy to Load onto the Car Trunk

The total width of the (8-inch) walker is 270mm, the depth is 750mm, and the height is 520mm when compactly folded. It is easy to store in the car trunk when going on drives.



Options Exclusive to OLive

**Walking Cane Holder (Small)** Suggested Retail Price (Non-Taxable) **¥5,000UP**

This can hold one walking cane. Please take the walking cane out of the holder before folding its handgrips backwards.

**Resistor** Suggested Retail Price (Non-Taxable) **¥9,000UP**

The speed can be adjusted depending on walking conditions for safety.


**Tipping-Less** **No Additional Fees**

This is option is offered by the manufacturer. (This part is not sold separately. Please ask for this option at the time of order.)

Standard Model (Fixed Backrest)



Seat Upholstery Color/Red



Standard Model (Fixed Backrest)



Seat Upholstery Color/Aqua Blue



New Concept - Compact Wheelchair.

**MV-2** Discretionary Price Caster x Rear Wheel [Inches] **4×6**

Front Seat Height **470 mm** Weight **6.2kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Rear Seat Height	Armrest Height	Backrest Height	Footrest Length	Total Width	Total Length	Total Height	Total Width (Folded)	Total Height (Folded)
<b>MV-2</b>	<b>470</b>	<b>350</b>	450	220	330	315	445	730	915	310	610

© Specifications are subject to change without prior notice for continual product improvement.



Seat Upholstery Color/Beige

- Ultra-lightweight! Weighs Only 6.2kg**
  - This can be used anywhere with its ultra-compact design. It is especially useful in public spaces and narrow entrances of bullet trains and buses. It is the perfect companion for outings.
  - All it takes is one touch to fold the wheelchair up. Its one-touch footbrake is a standard feature.
  - The highly practical casters rotate smoothly even on uneven terrain.

Safe and Comfortable. An All-Around Walking Aid for a Tireless Walk.

**MV-10** (Tax Included) **¥40,700** Caster x Rear Wheel [Inches] **4×4**

Front Seat Height **440 mm** Weight **3.5kg**

Unit: mm

Product Number	Front Seat Height	Seat Width	Rear Seat Height	Total Width	Total Length	Total Height	Total Width (Folded)	Total Height (Folded)
<b>MV-10</b>	<b>440</b>	<b>270</b>	405	390	615	780~925	240	860



© Specifications are subject to change without prior notice for continual product improvement.



Seat Upholstery Color/Broad Flowers

- Ultra Lightweight! Only 3.5kg for One-Handed Lifting**
  - Aluminum pipes are used for the main frame. It is painted with our company's exclusive surface coating for a high quality finish.
  - All it takes is one touch to fold the walking aid up for easy storage. It doubles as a very convenient chair for when the user is tired.
  - The bag for storing small items is a standard feature.

\* This is not a wheelchair. Please do not move it with a person seated on top.





4-Wheel • 6-Wheel Walker

A Safe and Reliable Item to  
Take Steps Forward.



SM-30 (4-Wheeler)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥61,600**

Weight  
12.3kg

Unit: mm

Product Number	① Inner Width of U-Shape	② Inner Depth of U-Shape	③ Total Width	④ Total Length	⑤ Total Height	⑥ Total Width (Folded)	⑦ Total Height (Folded)	⑧ Turning Radius
SM-30	330~430	425	550~650	725	860~1,200	260	1,075	755~860

◎ Specifications are subject to change without prior notice for continual product improvement.

\* ⑧ Turning radius...755mm when ① is 330mm, or 860mm when ① is 430mm



SM-35 (6-Wheeler)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥70,400**


Weight  
13.3kg

Unit: mm

Product Number	① Inner Width of U-Shape	② Inner Depth of U-Shape	③ Total Width	④ Total Length	⑤ Total Height	⑥ Total Width (Folded)	⑦ Total Height (Folded)	⑧ Turning Radius
SM-35	330~430	425	550~650	725	860~1,200	260	1,075	445~475

◎ Specifications are subject to change without prior notice for continual product improvement.

\* ⑧ Turning radius...445mm when ① is 330mm, or 475mm when ① is 430mm



Folded

This model can be converted to a 4-wheel model to use as a conventional walker.



SM-40 (6-Wheeler with Brakes)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥77,000**

Weight  
14.4kg

Unit: mm


Product Number	① Inner Width of U-Shape	② Inner Depth of U-Shape	③ Total Width	④ Total Length	⑤ Total Height	⑥ Total Width (Folded)	⑦ Total Height (Folded)	⑧ Turning Radius
SM-40	330~430	425	550~650	725	860~1,200	260	1,075	445~475

◎ Specifications are subject to change without prior notice for continual product improvement.


\* ⑧ Turning radius...445mm when ① is 330mm, or 475mm when ① is 430mm

Specifications

Folded




Turning Radius 445mm



6-Wheel Walker (SM-35 - 40)  
Inner Shape of U-Shape at 330mm

Turning Radius 640mm



Conventional Model  
(SM-10 - 4-Wheeler)  
Inner Width of U-Shape at 320mm

Feature. 1

Resistor

The resistor can be fixed onto the walker using screws. It can be adjusted according to the physical condition of the user for safe use under optimal settings.



Feature. 2

Width Adjustment

Adjust the width of the walker by raising or lowering the X-shaped section. This mechanism is 6-level adjustable for use by a wide range of users from shorter to taller users. (Inner Width of U-Shape 330~430mm)



Feature. 3

Height Adjustment

Similar to conventional walkers, the U-shape section is height adjustable to an optimal setting for the user's body size and physical condition.



Reciprocal Walking Frame/Standard Walking Frame

Anterior • Posterior Rubber Tips Model



CMS-90A

Folded

CMS-90A

(Reciprocal Walking Frame)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥18,000**


Weight  
2.6kg

CMS-90B

(Standard Walking Frame)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥16,800**

Weight  
2.6kg



Unit: mm

Product Number		Width Between Handgrips	Total Width	Total Length	Total Height	Diameter (Rubber Tips)	Color
CMS-90A	Reciprocal	460	550	470~500	730~890	19	Silver
CMS-90B	Standard	460	550	470~500	730~890	19	Silver

◎ Specifications are subject to change without prior notice for continual product improvement.

Anterior Swivel Caster Posterior Rubber Tip Model



CMS-91B

CMS-91A

(Reciprocal Walking Frame)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥21,600**


Weight  
2.8kg

CMS-91B

(Standard Walking Frame)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥20,400**

Weight  
2.8kg



Unit: mm

Product Number		Width Between Handgrips	Total Width	Total Length	Total Height	Diameter (Rubber Tips)	Color
CMS-91A	Reciprocal	460	550	510~540	800~960	19	Silver
CMS-91B	Standard	460	550	510~540	800~960	19	Silver

◎ Specifications are subject to change without prior notice for continual product improvement.

Anterior Swivel Caster Posterior Elastic Caster Model



CMS-92B

CMS-92A

(Reciprocal Walking Frame)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥24,000**


Weight  
2.9kg

CMS-92B

(Standard Walking Frame)

Manufacturer's Suggested Retail Price  
(Non-Taxable) **¥22,800**

Weight  
2.9kg



Unit: mm

Product Number		Width Between Handgrips	Total Width	Total Length	Total Height	Color
CMS-92A	Reciprocal	460	550	510~540	800~960	Silver
CMS-92B	Standard	460	550	510~540	800~960	Silver

◎ Specifications are subject to change without prior notice for continual product improvement.

Rubber Tips for Walking Canes

(Crutches, Quad Canes, Foldable Cane, Single Point Canes · Height Adjustable Walking Canes, Lofstrand Crutches, Seven Crutches)

φ16

φ19



A



B



C

These rubber tips are for crutches and quad canes.  
Two sizes for φ16 or φ19 tips are available.  
A...CMS-40S/M/L  
B...CMS-60S/L · CMS-70S/M/L  
C...CMS-80S/M/L

Rotating e-Pins

- Uses Rotating e-Pins
- Height Adjustable & Easy to Operate


ePin



Insert the pin inside the hole of the pipe and rotate to lock securely.



Walking Cane CMS Series



Crutches (2 per Set)			
Product Number	CMS-80L	CMS-80M	CMS-80S
Specification	Size L	Size M	Size S
Length (Total Length)	1,175~1,475 (13-Level)	1,075~1,375 (13-Level)	965~1,115 (7-Level)
Length (Underarm ~ Handgrip)	355~530 (8-Level)	305~455 (7-Level)	255~405 (7-Level)
Weight (Each)	1,000g	960g	940g
Color	Silver/Black		
Diameter (Rubber Tip)	19mm (Φ)		
Manufacturer's Suggested Retail Price	¥10,800 (Non-Taxable)		

© Specifications are subject to change without prior notice for continual product improvement.



Unit: mm					
	Lofstrand Crutches		Seven Crutches		
Product Number	CMS-60L	CMS-60S	CMS-70L	CMS-70M	CMS-70S
Specification	Size L	Size S	Size L	Size M	Size S
Length (Grip Height)	645~970	520~770	710~960	535~735	420~545
Length (Cuff Length)	200~300	180~260	200~300	180~260	160~220
Cuff Diameter	90mm(L)	70mm(M)	90mm(L)	70mm(M)	65mm(S)
Weight	630g	510g	595g	530g	480g
Color	Silver/Black		Silver/Black		
Diameter (Rubber Tip)	19(Φ)		19(Φ)		
Manufacturer's Suggested Retail Price	¥9,360 (Non-Taxable)		¥10,560 (Non-Taxable)		

© Specifications are subject to change without prior notice for continual product improvement.



Quad Cane			
Product Number	CMS-40L	CMS-40M	CMS-40S
Specification	Size L	Size M	Size S
Length	710~935 (10-Level)		
Base Size	315 (290) × 220	265 (250) × 190	220 (205) × 165
Weight	1,060g	1,020g	980g
Color	Silver/Black		
Diameter (Rubber Tip)	16mm (Φ)		
Manufacturer's Suggested Retail Price	¥7,680 (Non-Taxable)	¥7,440 (Non-Taxable)	¥7,200 (Non-Taxable)

© Specifications are subject to change without prior notice for continual product improvement.

Bathing Aids Coral Series

\* No demo units are available.



For barrier-free spaces...  
Please use  
a standard model

It features tight turns for  
easy use in narrow spaces.

SW-11 (Tax Included) ¥151,800

For spaces with uneven flooring like thresholds...  
Please use  
a deluxe model

For wider spaces, we recommend  
W Models with rear wheels fixed on the exterior.

SW-11W (Tax Included) ¥165,000

For narrower spaces, we recommend  
S Models with rear wheels fixed on the interior.

SW-11S (Tax Included) ¥165,000

Shower Carry

Toilet Carry

Weight 11.3kg

Weight 12.5kg

Weight 12.5kg

Weight 11.6kg

Weight 12.7kg

Weight 12.7kg

This can also be used for showering.

SW-21 (Tax Included) ¥158,400

This can also be used for showering.

SW-21W (Tax Included) ¥171,600

This can also be used for showering.

SW-21S (Tax Included) ¥171,600

The seat height of the toilet carry is adjustable.

Front Handgrip

For Deluxe Models Only

This aids users with unstable sitting positions to maintain proper sitting positions.

7-Inch Casters


For Deluxe Models Only

These casters are for rolling over thresholds. It is also equipped with a convenient tipping level and brakes to lock the wheels.

Adjustable Seat Height

For Toilet Carries Only

The height can be adjusted according to the height of the toilet seat.



Coral Series							
Main Body	Aluminum	Pad	Polyurethane Foam	Wheel	(Standard) All 4-Inch Wheels	Turning Radius	Standard: 475
Surface Finish	Powder Coating	Seat Pad	Width 440 × Depth 410		(Delux) Anterior 4-Inch Wheels		Deluxe: 725
Metal Part	Stainless	Backrest Pad	Width 400 × Height 200		Posterior 7-Inch Wheels		

Unit: mm



Shower Chair

Seat height is 8-level adjustable from 370 to 510mm at 20mm increments

Weight  
2.1kg



SC-11 (Tax Included) ¥10,890

Weight  
2.7kg



SC-12 (Tax Included) ¥14,520

Weight  
3.2kg



SC-13 (Tax Included) ¥15,730

Weight  
3.2kg



SC-21 (Tax Included) ¥13,310

Weight  
3.8kg



SC-22 (Tax Included) ¥16,940

Weight  
4.5kg



SC-23 (Tax Included) ¥18,150

\* Please refer to the user's manual available on Matsunaga Manfactory's website for the dimensions of each model.

Handgrips

Select from three types of handgrips: high handgrips · low handgrips · no handgrips.

Low Handgrips



1 per Side

The handgrips are installed at heights easy to grab for users with rounded backs.

High Handgrips



1 per Side

The high handgrips help maintain good posture. Easy to grip, the surface is treated with non-slip coating.

e-pin

Use the e-pins to securely lock into place to adjust the height.



Chair Legs  
4 Spots

Insert the pin inside the hole of the pipe and rotate to lock securely.

Shower Hook Included

The shower head can be hooked into place.



1 per Side

This frees up both hands when washing the body.

Foldable Shower Chair

Seat height is 6-level adjustable from 350 to 450mm at 20mm increments.

- The chair is safe to use, locking open when unfolded.
- The chair is easy to fold up.



SC-31 (Tax Included) ¥22,748

Specifications	SC-31
Total Width	495
Depth	450~550
Total Height	680~780
Seat Height	350~450
Dimensions (Folded)	170

Unit: mm

Bathtub Chair

The seat height is 4-level adjustable from 170 to 230mm at 20mm increments.

Detachable & Replacable Non-Slip Rubber Seat

In a Color Easy to Spot in Steam/Water

e-Pin

Uses e-Pins that Lock Securely to Adjust Height

Legs  
4 Spots

Suction Disc Model

Weight  
2.5kg



YC-1 (Tax Included) ¥15,730

Rubber Tips Model


Weight  
2.4kg



YC-2 (Tax Included) ¥14,520


Easily Detachable

1 per Side



Legs  
4 Spots


A very safe and secure design using rubber tips with double-layered suction discs.



YC-1	YC-2
Material/Frame: Stainless	
Top: Polypropylene (Antibacterial)	
Rubber Seat: Elastomer (Antibacterial)	

Handrail for Bathtubs

- The part which comes into contact with the bathtub is equipped with rubber sheets to prevent damages.
- Bathtub lids can be used with the handrail attached.
- The handrail comes with 3 extra boards to fill any gaps.



YG-1 (Tax Included) ¥22,990

Specifications	YG-1
Main Body Dimension	W 210 × D 225~285 × H 315~395
Height of Handrail to Get In/Out	160~240
Suitable Bathtub Rim Width	40~150
Handrail for Inside the Bathtub	Attachable on Left or Right Side
Frame	Stainless
Resin Parts	Polypropylene

Unit: mm



## NOTE

---



## Wheeled Stretcher & Stretcher





Aluminum Lightweight Wheeled Stretcher

Lightweight Wheeled Stretcher (High-Low Model)



Crank  
Handle  
Height  
Adjustable


Fixed  
Stretcher

Easy to operate for all users.

AL-UD-1 (4-Wheel Double Lock Model)

(Tax Included) ¥300,080

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 525~910mm
- Casters/φ150mm Double Lock
- Weight/40.2kg
- Method of Height Adjustment/Crank Handle
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm




Crank  
Handle  
Height  
Adjustable

Fixed  
Stretcher

AL-TUD-1 (Total Lock Model)

(Tax Included) ¥392,040

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 525~910mm
- Casters/φ150mm Total Lock U-Shaped Model
- Weight/42.8kg
- Method of Height Adjustment/Crank Handle
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm



Crank  
Handle  
Height  
Adjustable


Fixed  
Stretcher

\* Pediatric Model of AL-UD-1

AL-UD-1S (4-Wheel Double Lock Model - Pediatric)


(Tax Included) ¥371,470

- Material/Aluminum
- Dimension/Total Width 500 x Total Length 1,500 x Total Height 610~920mm
- Casters/φ150mm Double Lock
- Weight/34.0kg
- Method of Height Adjustment/Crank Handle
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 390 x Total Length 1,300 x Total Height 50mm



What are double lock casters • total lock casters?

Double Lock Casters



Standard  
φ150mm

Lowering the lock lever simultaneously locks swivel movement **A** and rotational movement **B** of the casters.

Total Lock U-Shaped Model



Standard  
φ150mm

The following 2 locking patterns are engaged by raising and lowering the lock lever.  
① Simultaneously locks movements **A** and **B** for all 4 wheels.  
② Locks movement **A** of only 1 caster.  
This stretcher is useful for transporting patients in straight lines over long distances, e.g. in hospital hallways.

Standard Feature  
Irrigator Stand Set



Lightweight Wheeled Stretcher (High-Low Model - Raising Backrest)



Crank  
Handle  
Height  
Adjustable

Fixed  
Stretcher


8-Level  
Backrest  
Adjustable  
Mechanism  
Included

A Lightweight and Reliable  
Wheeled Stretcher.

AL-UD-2 (4-Wheel Double Lock Model)

(Tax Included) ¥338,800

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 565~945mm
- Casters/φ150mm Double Lock
- Weight/41.2kg
- Method of Height Adjustment/Crank Handle, Raising Backrest with Pin Style Jack
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,760 x Total Height 50mm



Crank  
Handle  
Height  
Adjustable


Fixed  
Stretcher

8-Level  
Backrest  
Adjustable  
Mechanism  
Included

AL-TUD-2 (Total Lock Model)

(Tax Included) ¥430,760

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 565~945mm
- Casters/φ150mm Total Lock U-Shaped Model
- Weight/43.8kg
- Method of Height Adjustment/Crank Handle, Raising Backrest with Pin Style Jack
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,760 x Total Height 50mm



Lightweight Wheeled Stretcher (High-Low Model - Removable Stretcher)




Crank  
Handle  
Height  
Adjustable

Removable  
Stretcher

AL-UD-3 (4-Wheel Double Lock Model)

(Tax Included) ¥319,440

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 565~940mm
- Casters/φ150mm Double Lock
- Weight/43.6kg
- Method of Height Adjustment/Crank Handle, Stretcher Lever Lock
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm




Crank  
Handle  
Height  
Adjustable

Removable  
Stretcher

AL-TUD-3 (Total Lock Model)

(Tax Included) ¥411,400

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 565~940mm
- Casters/φ150mm Total Lock U-Shaped Model
- Weight/46.2kg
- Method of Height Adjustment/Crank Handle, Stretcher Lever Lock
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm





Aluminum Lightweight Wheeled Stretcher


Lightweight Wheeled Stretchers (High-Low Model - Raising Backrest - Removable Stretcher)



- Crank Handle Height Adjustable
- Removable Stretcher
- 8-Level Backrest Adjustable Mechanism Included

AL-UD-4 (4-Wheel Double Lock Model) (Tax Included) ¥358,160


- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 595~975mm
- Casters/φ150mm Double Lock
- Weight/44.6kg
- Method of Height Adjustment/Crank Handle, Stretcher Lever Lock
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,760 x Total Height 50mm



- Crank Handle Height Adjustable
- Removable Stretcher
- 8-Level Backrest Adjustable Mechanism Included

AL-TUD-4 (Total Lock Model) (Tax Included) ¥450,120

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 595~975mm
- Casters/φ150mm Total Lock U-Shaped Model
- Weight/47.2kg
- Method of Height Adjustment/Crank Handle, Stretcher Lever Lock
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,760 x Total Height 50mm




Lightweight Wheeled Stretcher (for Showers)



- Crank Handle Height Adjustable
- Fixed Stretcher

AL-SUD (4-Wheel Double Lock Model - for Showers) (Tax Included) ¥338,800

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 495~875mm
- Casters/φ125mm Double Lock
- Weight/37.6kg
- Method of Height Adjustment/Crank Handle
- Standard Features/Polyethylene Foam, Irrigator Stand
- Mattress Dimensions/Total Width 490 x Total Length 1,700 x Total Height 30mm




Other Aluminum General-Purpose Models



- Fixed Stretcher

AL-D (4-Wheel Double Lock Model) (Tax Included) ¥139,150


- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 790mm
- Casters/φ150mm Double Lock
- Weight/26.9kg
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm



- Fixed Stretcher
- 8-Level Backrest Adjustable Mechanism Included

AL-B (4-Wheel Double Lock Model - Backrest Adjustable Mechanism Included) (Tax Included) ¥179,080

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 825mm
- Casters/φ150mm Double Lock
- Weight/27.9kg
- Method of Backrest Adjustment/Raising Backrest with Pin Style Jack
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,760 x Total Height 50mm



Compact Wheeled Stretcher



- Fixed Stretcher
- A Compact Wheeled Stretcher Taking Up Minimal Storage Space

Foldable Stretcher (Tax Included) ¥206,910

- Material/Aluminum
- Dimension/Total Width 585 x Total Length 1,900 x Total Height 785mm
- Folded Dimension/Total Width 585 x Total Length 420 x Total Height 1,195mm
- Casters/φ150mm Double Lock
- Weight/25.9kg
- Standard Features/Mattress, Irrigator Stand




\* Folded

\* The lock uses a safety lock mechanism to prioritize safety. It can be easily folded.

What are double lock casters • total lock casters?

Double Lock Casters



Standard φ150mm

Lowering the lock lever simultaneously locks swivel movement **A** and rotational movement **B** of the casters.

Total Lock U-Shaped Model



Standard φ150mm


The following 2 locking patterns are engaged by raising and lowering the lock lever.  
① Simultaneously locks movements **A** and **B** for all 4 wheels.  
② Locks movement **A** of only 1 caster. This stretcher is useful for transporting patients in straight lines over long distances, e.g. in hospital hallways.

Standard Feature Irrigator Stand Set



Common Options for Wheeled Stretchers \* Unavailable for SUD • MRI Series (Except Fixation Strap for Stretchers)

Oxygen Tank Holder



For 500cc Oxygen Tanks  
\* The photo is a model for AL-UD • TUD. Oxygen tank holders for attachment onto other types of wheeled stretchers require alterations to the main body.

Basket



Total Width 450 x Total Length 700 x Total Height 100 mm  
Surface coated

Fixation Strap for Wheeled Stretchers



L 1,900 x W 70 mm



MRI-Compatible Wheeled Stretcher

3 Tesla MRI-Compatible Wheeled Stretcher



Seat Color/B-6

Crank  
Handle  
Height  
Adjustable


Fixed  
Stretcher

Non-magnetic aluminum material is used in the main frame.  
It is safe to use in 3 Tesla MRI rooms.

MR-UD-3T (4-Wheel Double Lock Model)

(Tax Included) ¥453,200

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 495-860mm
- Casters/φ125mm Double Lock, Stainless
- Weight/41.3kg
- Method of Height Adjustment/Crank Handle
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm



Please note that alterations are unavailable for this model.

Ultra-High Spec 3 Tesla MRI-Compatible Wheeled Stretcher



Seat Color/B-38

Crank  
Handle  
Height  
Adjustable


Fixed  
Stretcher

The frame, bolt, nuts and other detailed parts are made of non-magnetic material. This minimizes the risk of disturbances to MRI images even within 3 Tesla MRI rooms.

MR-UD-3TMG (4-Wheel Double Lock Model)

(Tax Included) ¥968,000

- Material/Aluminum
- Dimension/Total Width 590 x Total Length 1,900 x Total Height 525-900mm
- Casters/φ150mm Double Lock, Made of Resin
- Weight/39.1kg
- Method of Height Adjustment/Crank Handle
- Standard Features/Mattress, Irrigator Stand
- Mattress Dimensions/Total Width 500 x Total Length 1,690 x Total Height 50mm



Please note that alterations are unavailable for this model.



MR-3T







MR-3TMG

\* Wheelchairs are also available.

Stretcher

Pole Stretcher



Storage Bag for 2-Fold Stretchers

Discretionary Price

Storage Bag for 4-Fold Stretchers


Discretionary Price

\* The stretcher in the photo is the 4-fold model.

2-Fold Stretcher with Legs (Fixed Handles) Made of Steel (Tax Included) ¥15,125

Made of Aluminium (Tax Included) ¥26,620

Specifications	Steel	Aluminum
Dimensions	Total Width 540 x Total Length 2,040mm	
Weight	7.0kg	5.4kg




The most economical option out of our selection of standard stretchers.

2-Fold Stretcher with Legs (Retractable Handles) Made of Steel (Tax Included) ¥22,385

Made of Aluminium (Tax Included) ¥29,040

Specifications	Steel	Aluminum
Dimensions	Total Width 540 x Total Length 2,270mm Rolled Up Total Length 1,790mm	
Weight	8.5kg	6.4kg




The handles are retractable for use in narrow spaces and inside elevators.

4-Fold Stretcher with Legs (Fixed Handles) Made of Steel (Tax Included) ¥25,410

Made of Aluminium (Tax Included) ¥30,250

Specifications	Steel	Aluminum
Dimensions	Total Width 580 x Total Length 2,040mm 4-Fold Total Length 1,060mm	
Weight	8.1kg	6.5kg




The stretcher can be compactly folded into 4.

4-Fold Stretcher with Legs (Retractable Handles) Made of Steel (Tax Included) ¥27,830

Made of Aluminium (Tax Included) ¥32,670

Specifications	Steel	Aluminum
Dimensions	Total Width 580 x Total Length 2,270mm Rolled Up Total Length 1,790mm 4-Fold Total Length 920mm	
Weight	9.6kg	7.5kg



The handles are retractable and easy to grab.

Options

● **Stretchers**

- 2-Fold Stretcher with Legs & Fixed Handles
- 2-Fold Stretcher with Legs & Retractable Handles
- 4-Fold Stretcher with Legs & Fixed Handles
- 4-Fold Stretcher with Legs & Retractable Handles

● **Seat**

Tarpaulin

+¥4,400 (Tax Included)

● **Seat**

Made of Vinyl Coating

+¥2,200 (Tax Included)

● **Casters**

They can be attached to one side of the stretcher for smooth operation.

+¥6,600 (Tax Included)

● **Seat**

Dark Green Waterproof Canvas

+¥2,200 (Tax Included)


Wheeled  
Stretcher  
Stretcher

Stretcher

Wheeled  
Stretcher  
Stretcher

MRI-Compatible Wheeled Stretcher





## A Model Stretcher 1

Aluminum (Tax Included) **¥94,380**


Specification	Aluminum
Dimension	Total Width 515 x Total Length 1,900 x Total Height 215mm    Folded Total Length 960mm
Weight	8.0kg



\* Folded

Casters are attached on one side of the stretcher to be carried by 1 person.

It can be transported by one person. The post crutch and casters are retractable with one-touch. The stretcher can be compactly folded for storage.



## A Model Stretcher 2

Aluminum (Tax Included) **¥78,650**


Specification	Aluminum
Dimension	Total Width 490 x Total Length 1,900 x Total Height 60mm    Folded Total Length 960mm
Weight	5.0kg



\* Folded

It can be folded in half for storage.

It is lightweight and can be compactly stored. It can be taken on outdoor activities.



\* Chair Mode


## Chair-Type Stretcher

Aluminum (Tax Included) **¥118,580**

Specification	Aluminum
Dimension	Total Width 550 x Total Length 1,925 x Total Height 255mm Chair Mode : Total Width 550 x Total Length 650 x Total Height 1,395mm
Weight	About 8.5kg

A cart/stretcher which can be transported by one person.

One person can easily carry this stretcher.



## Assist Stretcher

Aluminum (Tax Included) **¥166,980**

Specification	Aluminum		
	Flat Stretcher Mode	Chair Mode	Folded
Total Length	1,865mm	660mm	985mm
Total Width	520mm		505mm
Total Height	195mm	1,360mm	140mm
Seat to Floor Height	—	520mm	—
Backrest Height	850mm		—
Weight	9.5kg		
Maximum Load Capacity	120kg		—

The Ultimate Multi-Purpose Stretcher!!



# Full of Helpful Information!! Matsunaga's Social Media!

## Matsunaga's Website



This is where customers can access our list of parts and electronic catalogues. Customers can download instruction manuals and make inquiries here.

## YouTube



We release a wide range of videos about our wheelchairs, from introductory videos of new products to step-by-step tutorials on how to adjust backrests.

## Instagram



New media contents by Matsunaga Manufactory! We regularly post information about our products and exhibitions.

## Facebook



We regularly post information about Matsunaga. Customers can easily check out our YouTube videos here.

## Twitter



This is where we tweet about wheelchair-related content. We might even post behind-the-scene photos of our factories here!

Even if the QR codes cannot be successfully read, it is easy to find us by typing "Matsunaga Manufactory" into the search box of each platform. Please check us out whenever you have some free time!



We are also expanding our MP series - our brand for active users. Please check it out.

## Website for MP Series



## YouTube



## Facebook



## Instagram





Product Name	Page
A	
Aerial Series	38~40
AL-B	83
AL-D	83
AL-SUD	82
AL-TUD—1	80
AL-TUD—2	81
AL-TUD—3	81
AL-TUD—4	82
AL-UD—1	80
AL-UD—1S	80
AL-UD—2	81
AL-UD—3	81
AL-UD—4	82
A Model Stretcher 1	86
A Model Stretcher 2	86
AR-101	28
AR-101BT	66
AR-111	28
AR-111BT	66
AR-201B	29
AR-201BT	66
AR-211B	29
AR-271B	29
AR-280	31
AR-301	30
AR-311	30
AR-371	30
AR-380	31
AR-501	32
AR-511B	32
AR-511BT	66
AR-601	32
AR-901	33, 34
AR-911S	35
Assist Stretcher	86
C	
Chair-Type Stretcher	86
CM-251HB	26
CM-261HB	26
CM-251HB-SE	26
CM-261HB-SE	27
CM-501HB	54
CM-541HB	54
CMS-40L/M/S	74

Product Name	Page
CMS-60L/S	74
CMS-70L/M/S	74
CMS-80L/M/S	74
CMS-90A/90B	73
CMS-91A/91B	73
CMS-92A/92B	73
Coral Series	75~77
D	
DM-81	24
E	
Elise	57
Escort	56
F	
Foldable Stretcher	83
FR-11R	56
FR-31TR	57
G	
Grace Core Series	6~9
GRC-11B	6
GRC-21B	6
GRC-21BD	6
GRC-31B	7
GRC-41B	7
GRC-41BD	7
H	
Hop Step Series	72
L	
Leveret (MKD-M/L)	64
LUUK (-01/-02)	60
M	
MB-PONY	62
MH-CR3D	45
MH-GR/GR-SE/GRL/GRL-SE	46, 47
MH-SR/SR-SE/SRL/SRL-SE	44
MH-RD3D	48
MH-VR	49
MH-VR-SE	49
Mignon (Tilting Stroller)	63
MKD-M/L (Leveret)	64
morph	4, 42
MR-3T	67
MR-3TMG	67
MR-IG-3TMG	67
MR-UD-3T	84
MR-UD-3TMG	84


Product Name	Page
MV-2	71
MV-10	71
MV-100	70
MW-14HB	55
MW-SL11B	38
MW-SL21B	38
MW-SL31B	39
MW-SL41B	39
MW-SL5B	40
MW-SL5BT	66
MW-SL6B	40
My Tilt - Compact 3D	45
My Tilt - Grand 3D	46, 47
My Tilt - Mini 3D	44
My Tilt - Rigid 3D	48
My Tilt - Value	49
MY-1	36
MY-2	36
N	
NEW MB-1	62
NEXT CORE Series	11~22
NEXT-10B HB/10B	17
NEXT-11B HB/11B	16
NEXT-12B HB/12B	17
NEXT-20B HB/20B	17
NEXT-21B HB/21B/21B HBD	16
NEXT-22B HB/22B	17
NEXT-31B HB/31B	15
NEXT-41B HB/41B/41B HBD	15
NEXT-50B HB/50B	13
NEXT-51B HB/51B	12
NEXT-52B HB	14
NEXT-60B HB/60B	13
NEXT-61B HB/61B/61B HBD	12
NEXT-62B HB	14
NEXT-70B HB/70B	18
NEXT-71B HB/71B	18
NEXT-80B HB/80B	18
NEXT-81B HB/81B	18
O	
OASIS Series	50~53
OH-1	68
OLive	70
One-Arm Drive Wheelchair	68
OS-12TR	52

Product Name	Page
OS-11TRSP	51
OS-12TRSP	51
S	
SC-11	76
SC-12	76
SC-13	76
SC-21	76
SC-22	76
SC-23	76
SC-31	76
SM-30	72
SM-35	72
SM-40	72
Stretchers	79
-Aluminum Wheeled	80
-MRI-Compatible	84
-Pole series	85
-Others	86
SW-11/11W/11S	75
SW-21/21W/21S	75
T	
TT-Jr	64
Y	
YC-1	77
YC-2	77
YG-1	77





Company Name	Matsunaga Manufactory Co., Ltd.
President	Noriyuki Matsunaga
Location	484 Ohba, Yoro-cho, Yoro-gun, Gifu, Japan
Capital	¥50,000,000 JPY
Gross Yearly Sales	¥6,000,000,000 JPY (Based on May 2020 Financial Results)
Number of Employees	179
Plot Size of Head Factory	18,876m <sup>2</sup>
Floor Space of Head Factory	17,820m <sup>2</sup>



\* This catalogue is the edited version of Japanese catalogue Vol.19-1 and Vol.20-1.





Head Office



Head Office Building



Shanghai Office Building



Shanghai Factory




Leather Processing Machine



3D Bending Machine



3D Router



Welding Line



Delivery Line





MATSUNAGA

## Matsunaga Manufactory Co., Ltd.

### ■ Head Office

484 Ohba, Yoro-cho, Yoro-gun, Gifu, 503-1272, Japan  
tel. +81 584-35-1180 (Main) fax. +81-584-35-1270

### ■ Tokyo Sales Office

1-19-12 Minami-Shinozakimachi, Edogawa, Tokyo, 133-0065, Japan  
tel. +81-3-5636-1180

U <https://www.matsunaga-w.co.jp/global/>

✉ [info@matsunaga-w.co.jp](mailto:info@matsunaga-w.co.jp)

- The catalog price is only valid for Japan domestic market from May 2022.  
Products listed as non-taxable are exempt from consumption tax, etc.
- Specifications and prices are subject to change without prior notice for product development and improvement.
- Please note that the colors in product pictures and color samples may vary slightly from the actual colors of the products due to printing.