

* No demo units are available.

For barrier-free spaces...

Please use a standard model

It features tight turns for easy use in narrow spaces.

For spaces with uneven flooring like thresholds...

Please use a deluxe model

For wider spaces, we recommend W Models with rear wheels fixed on the exterior.

For narrower spaces, we recommend S Models with rear wheels fixed on the interior.

Shower Carry



Weight 11.3kg

SW-11 (Tax Included) ¥151,800



Weight 12.5kg

SW-11W (Tax Included) ¥165,000



Weight 12.5kg

SW-11S (Tax Included) ¥165,000

Toilet Carry



Weight 11.6kg

This can also be used for showering.
SW-21 (Tax Included) ¥158,400



Weight 12.7kg

This can also be used for showering.
SW-21W (Tax Included) ¥171,600



Weight 12.7kg

This can also be used for showering.
SW-21S (Tax Included) ¥171,600

The seat height of the toilet carry is adjustable.

Front Handgrip

For Deluxe Models Only



This aids users with unstable sitting positions to maintain proper sitting positions.

7-Inch Casters

For Deluxe Models Only



These casters are for rolling over thresholds. It is also equipped with a convenient tipping level and brakes to lock the wheels.



Adjustable Seat Height

For Toilet Carries Only





The height can be adjusted according to the height of the toilet seat.

■ SW-11



■ SW-21W/S



● Coral Series

| | | | | | | | |
|----------------|----------------|--------------|------------------------|-------|--------------------------------|----------------|---------------|
| Main Body | Aluminum | Pad | Polyurethane Foam | Wheel | (Standard) All 4-Inch Wheels | Turning Radius | Standard: 475 |
| Surface Finish | Powder Coating | Seat Pad | Width 440 × Depth 410 | | (Delux) Anterior 4-Inch Wheels | | Deluxe: 725 |
| Metal Part | Stainless | Backrest Pad | Width 400 × Height 200 | | Posterior 7-Inch Wheels | | |

Unit: mm

Shower Chair

Seat height is 8-level adjustable from 370 to 510mm at 20mm increments

Weight 2.1kg



SC-11 (Tax Included) ¥10,890

Weight 2.7kg



SC-12 (Tax Included) ¥14,520

Weight 3.2kg



SC-13 (Tax Included) ¥15,730

Weight 3.2kg



SC-21 (Tax Included) ¥13,310

Weight 3.8kg



SC-22 (Tax Included) ¥16,940

Weight 4.5kg



SC-23 (Tax Included) ¥18,150

* Please refer to the user's manual available on Matsunaga Manfactory's website for the dimensions of each model.

Handgrips

Select from three types of handgrips: high handgrips · low handgrips · no handgrips.

Low Handgrips



1 per Side

The handgrips are installed at heights easy to grab for users with rounded backs.

High Handgrips



1 per Side

The high handgrips help maintain good posture. Easy to grip, the surface is treated with non-slip coating.

e-pin

Use the e-pins to securely lock into place to adjust the height.



Chair Legs 4 Spots

Insert the pin inside the hole of the pipe and rotate to lock securely.

Shower Hook Included

The shower head can be hooked into place.



1 per Side

This frees up both hands when washing the body.

Foldable Shower Chair

Seat height is 6-level adjustable from 350 to 450mm at 20mm increments.

- The chair is safe to use, locking open when unfolded.
- The chair is easy to fold up.



Weight 4.1kg

Equipped with a Handle

In a Color Easy to Spot in Steam

Flip-Up Armrest They are convenient for easy transfers.

Equipped with anti-bacterial & skin-friendly soft cushions They can be removed and washed.

e-Pin e-Pins to lock securely into place

Chair Legs 4 Spots

SC-31 (Tax Included) ¥22,748

| Specifications | SC-31 |
|---------------------|---------|
| Total Width | 495 |
| Depth | 450~550 |
| Total Height | 680~780 |
| Seat Height | 350~450 |
| Dimensions (Folded) | 170 |

Bathtub Chair

The seat height is 4-level adjustable from 170 to 230mm at 20mm increments.

Detachable & Replacable Non-Slip Rubber Seat

In a Color Easy to Spot in Steam/Water

e-Pin Uses e-Pins that Lock Securely to Adjust Height

Legs 4 Spots

Suction Disc Model

Weight 2.5kg



YC-1 (Tax Included) ¥15,730

Rubber Tips Model

Weight 2.4kg




YC-2 (Tax Included) ¥14,520

Easily Detachable

1 per Side




Legs 4 Spots A very safe and secure design using rubber tips with double-layered suction discs.



| YC-1 | YC-2 |
|--|------|
| Material/Frame: Stainless | |
| Top: Polypropylene (Antibacterial) | |
| Rubber Seat: Elastomer (Antibacterial) | |

Handrail for Bathtubs

- The part which comes into contact with the bathtub is equipped with rubber sheets to prevent damages.
- Bathtub lids can be used with the handrail attached.
- The handrail comes with 3 extra boards to fill any gaps.



Handrail to Get In/Out The handgrip form is easy to grab from any direction. Height adjustable.

In a Color Easy to Spot in Steam

Easy to grab Non-slip coating

Handrail for Inside the Bathtub The handgrip form is comfortable to grab while taking a bath. It is attachable on the left or right side.

Weight 5.2kg



YG-1 (Tax Included) ¥22,990

| Specifications | YG-1 |
|----------------------------------|----------------------------------|
| Main Body Dimension | W 210 × D 225~285 × H 315~395 |
| Height of Handrail to Get In/Out | 160~240 |
| Suitable Bathtub Rim Width | 40~150 |
| Handrail for Inside the Bathtub | Attachable on Left or Right Side |
| Frame | Stainless |
| Resin Parts | Polypropylene |