


Crutches



The handgrip height is easy to adjust.



Insert the pin inside the hole of the pipe and rotate to lock securely.

CMS-80L CMS-80M CMS-80S Unit: mm

| Crutches (2 per Set) | | | |
|---------------------------------------|------------------------|------------------------|---------------------|
| Product Number | CMS-80L | CMS-80M | CMS-80S |
| Specification | Size L | Size M | Size S |
| Length (Total Length) | 1,175~1,475 (13-Level) | 1,075~1,375 (13-Level) | 965~1,115 (7-Level) |
| Length (Underarm ~ Handgrip) | 355~530 (8-Level) | 305~455 (7-Level) | 255~405 (7-Level) |
| Weight (Each) | 1,000g | 960g | 940g |
| Color | Silver/Black | | |
| Diameter (Rubber Tip) | 19mm (∅) | | |
| Manufacturer's Suggested Retail Price | ¥10,800 (Non-Taxable) | | |

© Specifications are subject to change without prior notice for continual product improvement.

Quad Cane



Strap (Cord) Included

Height Adjustable



Insert the pin inside the hole of the pipe and rotate to lock securely.

CMS-40L CMS-40M CMS-40S Unit: mm

| Quad Cane | | | |
|---------------------------------------|----------------------|----------------------|----------------------|
| Product Number | CMS-40L | CMS-40M | CMS-40S |
| Specification | Size L | Size M | Size S |
| Length | 710~935 (10-Level) | | |
| Base Size | 315 (290) × 220 | 265 (250) × 190 | 220 (205) × 165 |
| Weight | 1,060g | 1,020g | 980g |
| Color | Silver/Black | | |
| Diameter (Rubber Tip) | 16mm (∅) | | |
| Manufacturer's Suggested Retail Price | ¥7,680 (Non-Taxable) | ¥7,440 (Non-Taxable) | ¥7,200 (Non-Taxable) |

© Specifications are subject to change without prior notice for continual product improvement.

Lofstrand Crutches



Seven Crutches



Unit: mm

| Lofstrand Crutches | | Seven Crutches | | | |
|---------------------------------------|----------------------|----------------|-----------------------|----------|----------|
| Product Number | CMS-60L | CMS-60S | CMS-70L | CMS-70M | CMS-70S |
| Specification | Size L | Size S | Size L | Size M | Size S |
| Length (Grip Height) | 645~970 | 520~770 | 710~960 | 535~735 | 420~545 |
| Length (Cuff Length) | 200~300 | 180~260 | 200~300 | 180~260 | 160~220 |
| Cuff Diameter | 90mm (L) | 70mm (M) | 90mm (L) | 70mm (M) | 65mm (S) |
| Weight | 630g | 510g | 595g | 530g | 480g |
| Color | Silver/Black | | Silver/Black | | |
| Diameter (Rubber Tip) | 19 (∅) | | 19 (∅) | | |
| Manufacturer's Suggested Retail Price | ¥9,360 (Non-Taxable) | | ¥10,560 (Non-Taxable) | | |

© Specifications are subject to change without prior notice for continual product improvement.